VARIAN

Varian, Inc.
2700 Mitchell Drive
Walnut Creek, CA 94598-1675/usa

ProStar 430 AutoSampler

Operation Manual

VARIAN

Quality Systems At Varian, Inc.

The ISO 9000 series standards were created in Geneva in 1987 to cut through a morass of conflicting quality definitions. These standards define a model for quality assurance systems in product design, development, manufacturing, installation, service, and customer support. They are now the worldwide quality assurance benchmark used to gauge the strength of a company's commitment to quality, and the value of its quality systems.

Various organizations around the world, such as the British Standards Institution (BSI), provide certified, objective auditors to scrutinize quality procedures, product development, manufacturing processes, and customer satisfaction programs. No company can claim ISO 9000 series registration unless it receives a stamp of approval from the demanding quality assessors of BSI or similar accredited examining body. ISO 9000 series registration constitutes an objective third-party report to determine the level of a supplier's commitment to quality.

In 1992, Varian, Inc., Analytical Instruments became registered to the most comprehensive of the ISO 9000 series standards — ISO 9001. ISO 9001 registration means that every stage of our quality system, including product development, manufacturing, final test, shipping, and parts and supplies has been rigorously examined against the most exacting set of internationally recognized standards. It means we live up to a standard of quality that you can count on today, and into the future. Our Quality System has received ISO 9001 certification number FM21797.

The quality systems that earned us ISO 9001 registration have direct benefits for our customers:

- ◆ We can speed instruments to you faster than ever before. Emergency orders can be processed even faster.
- ◆ We fill your orders promptly and completely.
- ◆ We have implemented a system of continuous feedback from our customers — we are aware of your needs today and tomorrow.
- ◆ We have improved your productivity by cutting systems failure rates in half and speeding service response time.
- ◆ We have embedded continuous improvement into the fabric of our organization so that we can achieve even higher levels of quality in the future.
- ◆ We are embedding GLP requirements into our products and services to help you meet your regulatory compliance requirements.

ISO 9001 registration is not enough. For us, quality is defined by our customers. We are not satisfied unless you are satisfied. We are striving to understand customer needs, using independent surveys, user groups, customer advisory boards, and our "Hallmark of Quality" response program, in addition to individual face-to-face customer contact. Our products and our processes are configured to meet those needs.

We know that you are seeking more than the most advanced processes and top-notch applications expertise. You want to join forces with a partner committed to delivering world-class quality, reliability, and value — on time, every time.

Our overriding aim is to be that partner.

Varian, Inc. Analytical Instrument Warranty

Hardware Products

All analytical instruments sold by Varian, Inc. are warranted to be free from defects in material and workmanship for the periods specified and in accordance with the terms on the face of Varian's quotation or as otherwise agreed upon in writing between Varian and the Customer. The warranty period begins on the date of **shipment** from Varian to the original Customer. However, where installation is paid for by the Customer or included in the purchase price, the warranty period begins upon completion of installation. If the Customer schedules **installation** to start later than 30 days after delivery or if such delay is caused through the Customer's inability to provide adequate facilities or utilities or through failure to comply with Varian's reasonable pre-installation instructions or through other omissions by Customer, then the warranty period starts on the 31st day from date of shipment. Moreover Varian will charge the Customer for labor and other expenses involved in making multiple or follow-up installation service calls.

Software Products

Where software is provided within the frame of a license agreement concluded between the Customer and Varian, any warranty shall be strictly in accordance with the terms of such agreement.

In the absence of a license agreement and unless an alternate warranty period is agreed upon in writing between Varian and the Customer, the warranty period is as specified on the face of Varian's quotation. Varian warrants such software products, if used with and properly installed on Varian hardware or other hardware as specified by Varian to perform as described in the accompanying Operator's Manual and to be substantially free of those defects which cause failure to execute respective programming instructions; however, Varian does not warrant uninterrupted or error-free operation.

Remedies

The sole and exclusive remedy under hardware warranty shall be **repair** of instrument malfunctions which in Varian's opinion are due or traceable to defects in original materials or workmanship or, at Varian's option, **replacement** of the respective defective parts, provided that Varian may as an alternative elect to **refund** an equitable portion of the purchase price of the instrument or accessory.

Repair or replacement under warranty does not extend the original warranty period.

Repair or replacement under warranty claims shall be made in Varian's sole discretion either by sending a Customer Support Representative to the site or by authorizing the Customer to return the defective accessory or instrument to Varian or to send it to a designated service facility. The Customer shall be responsible for loss or damage in transit and shall prepay shipping cost. Varian will return the accessory or instrument to the Customer prepaid and insured. Claims for loss or damage in transit shall be filed by the Customer. To correct software operation anomalies, Varian will issue software revisions where such revisions exist and where, in Varian's opinion, this is the most efficient remedy.

Limitation of Warranty

This **warranty does not cover** software supplied by the Customer, equipment and software warranted by another manufacturer or replacement of expendable items and those of limited life, such as but not limited to: Filters, glassware, instrument status lamps, source lamps, septa, columns, fuses, chart paper and ink, nebulizers, flow cells, pistons, seals, fittings, valves, burners, sample tubes, probe inserts, print heads, glass lined tubing, pipe and tube fittings, variable temperature dewars, transfer lines, flexible discs, magnetic tape cassettes, electron multipliers, filaments, vacuum gaskets, seats and all parts exposed to samples and mobile phases.

This **warranty shall be void** in the event of accident, abuse, alteration, misuse, neglect, breakage, improper operation or maintenance, unauthorized or improper modifications or tampering, use in an unsuitable physical environment, use with a marginal power supply or use with other inadequate facilities or utilities. Reasonable care must be used to avoid hazards.

This warranty is expressly in lieu of and excludes all other express or implied warranties, including but not limited to warranties of merchantability and of fitness for particular purpose, use or application, and all other obligations or liabilities on the part of Varian, unless such other warranties, obligations or liabilities are expressly agreed to in writing by Varian.

Limitation of Remedies and Liability

The remedies provided herein are the sole and exclusive remedies of the Customer. In no case will Varian be liable for incidental or consequential damages, loss of use, loss of production or any other loss incurred.

Safety Information

Operating Instructions

This instruction manual is provided to help you establish operating conditions which will permit safe and efficient use of your equipment. Special considerations and precautions are also described in the manual, which appear in the form of **NOTES**, **CAUTIONS**, and **WARNINGS** as described below. It is important that you operate your equipment in accordance with this instruction manual and any additional information which may be provided by Varian. Address any questions regarding the safe and proper use of your equipment to your local Varian office.

NOTE

Information to aid you in obtaining optimal performance from your instrument.

CAUTION

Alerts you to situations that may cause moderate injury and/or equipment damage, and how to avoid these situations.

WARNING

Alerts you to potentially hazardous situations that could result in serious injury, and how to avoid these situations.

Warning Symbol

Warning Description

**WARNING:
SHOCK HAZARD**

Hazardous voltages are present inside instrument. Disconnect from main power before removing screw-attached panels.

**WARNING:
CHEMICAL HAZARD**

Hazardous chemicals may be present. Avoid contact, especially when replenishing reservoirs. Use proper eye and skin protection.

**WARNING:
BURN HAZARD**

Very hot or cryogenically cold surfaces may be exposed. Use proper skin protection.

**WARNING:
EYE HAZARD**

Eye damage could occur either from flying particles, chemicals, or UV radiation. Use proper eye and face protection.

**WARNING:
FIRE HAZARD**

The potential for fire may be present. Follow manual instructions for safe operation.

**WARNING:
EXPLOSION HAZARD**

The potential for explosion may exist because of type of gas or liquid used.

**WARNING:
RADIATION SOURCE**

Ionizing radiation source is present. Follow manual instructions for safe operation.

**WARNING:
MOVING PARTS**

Keep hands and fingers away.

General Safety Precautions

Follow these safety practices to ensure safe equipment operation.

- Perform periodic leak checks on all supply lines and pneumatic plumbing.
- Do not allow gas lines to become kinked or punctured. Place lines away from foot traffic and extreme heat or cold.
- Store organic solvents in fireproof, vented and clearly labeled cabinets so they are easily identified as toxic and/or flammable materials.
- Do not accumulate waste solvents. Dispose of such materials through a regulated disposal program and not through municipal sewage lines.

NOTICE: This instrument has been tested per applicable requirements of EMC Directive as required to carry the European Union CE Mark. As such, this equipment may be susceptible to radiation/interference levels or frequencies which are not within the tested limits.

WARNING

This instrument is designed for chromatographic analysis of appropriately prepared samples. It must be operated using appropriate gases and/or solvents and within specified maximum ranges for pressure, flows, and temperatures as described in this manual. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

WARNING

It is the responsibility of the Customer to inform Varian Customer Support Representatives if the instrument has been used for the analysis of hazardous biological, radioactive, or toxic samples, prior to any instrument service being performed or when an instrument is being returned to the Service Center for repair.

Electrical Hazards

- Disconnect the instrument from all power sources before removing protective panels to avoid exposure to potentially dangerous voltages.
- When it is necessary to use a non-original power cord plug, make sure the replacement cord adheres to the color coding and polarity described in the manual and all local building safety codes.
- Replace blown fuses with fuses of the size and rating stipulated on the fuse panel or in the manual.
- Replace faulty or frayed power cords immediately with the same type and rating.
- Make sure that voltage sources and line voltage match the value for which the instrument is wired.

Compressed Gas Cylinders

- Store and handle compressed gases carefully and in strict adherence to safety codes.
- Secure cylinders to an immovable structure or wall.
- Store and move cylinders in an upright, vertical position. Before transport, remove regulators and install cylinder cap.
- Store cylinders in a well-ventilated area away from heat, direct sunshine, freezing temperatures, and ignition sources.
- Mark cylinders clearly so there is no doubt as to their contents.
- Use only approved regulators and connections.
- Use only connector tubing that is chromatographically clean (Varian Part Number 03-918326-00) and has a pressure rating significantly greater than the highest outlet pressure from the regulator.

GC Safety Practices

Exhaust System

No special exhaust ducting is necessary for GC detectors installed in a well-ventilated room except when the detectors are used to test hazardous chemicals. If you do install ducting:

- Use only fireproof ducting.
- Install a blower at the duct outlet.
- Locate duct intakes such that their vibration or air movement does not effect detector operation.
- Check periodically for proper operation of the duct.
- Ensure proper ventilation in lab area.

Radioactive Source Detectors

- Read carefully and comply with all NOTES, CAUTIONS, and WARNINGS in the Ni⁶³ ECD manual.
- Perform the tests for removable radioactive contamination described in the Ni⁶³ ECD manual.
- Comply with leak test schedules and procedures.

Burn Hazard

Heated or cryogenically cooled zones of gas chromatographs can remain hot or cold for a considerable time after instrument power is turned off. To prevent painful burns, ensure that all heated or cooled areas have returned to room temperature or wear adequate hand protection before you touch potentially hot or cold surfaces.

LC Safety Practices

High Pressure Hazard

- If a line ruptures, a relief device opens, or a valve opens accidentally under pressure, potentially hazardous high liquid pressures can be generated by the pump causing a high velocity stream of volatile and/or toxic liquids.
- Wear face protection when you inject samples or perform routine maintenance.
- Never open a solvent line or valve under pressure. Stop the pump first and let the pressure drop to zero.
- Use shatter-proof reservoirs capable of operating at 50-60 psi.
- Keep the reservoir enclosure closed when the reservoir is under pressure.
- Read and adhere to all NOTES, CAUTIONS, and WARNINGS in the manual.

Flash Chromatography

The operator should be familiar with the physico-chemical properties of the components of the mobile phase.

Keep solvents from direct contact with the polyurethane supply tubing as certain solvents will cause weakening and leaks with possible bursting.

All components of the system should be connected to a common power supply and common ground. This ground must be a true ground rather than a floating ground.

Non-polar solvents can develop a static charge when pumped through the system. All vessels that contain mobile phase (including tubing and collection vessels) must be grounded to dissipate static electricity.

Employ static measuring and static discharge devices (e.g., air ionizers) to safeguard against the buildup of static electricity.

Ultraviolet Radiation

Liquid chromatograph detectors that use an ultraviolet light source have shielding to prevent radiation exposure to personnel.

For continued protection:

- Ensure that protective lamp covers of variable and fixed wavelength detectors are in place during operation.
- Do not look directly into detector fluid cells or at the UV light source. When inspecting the light source or fluid cell, always use protective eye covering such as borosilicate glass or polystyrene.

The following is a Federal Communications Commission advisory: This equipment has been tested and found to comply with the limits of a Class A computing device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Spare Parts Availability

It is the policy of Varian to provide operational spare parts for any instrument and major accessory for a period of five (5) years after shipment of the final production run of that instrument. Spare parts will be available after this five (5) year period but on an *as available* basis. Operational spare parts are defined as those individual electrical or mechanical parts that are susceptible to failure during their normal operation. Examples include relays, lamps, temperature probes, detector elements, motors, etc. Sheet metal parts, structural members or assemblies and castings, printed circuit boards, and functional modules are normally capable of being rebuilt to like-new condition throughout their useful life and therefore will be supplied only on an *as available* basis after the final production run of the instrument.

Service Availability

Varian provides a variety of services to support its customers after warranty expiration. Repair service can be provided by attractively priced service contracts or on a time and material basis. Technical support and training can be provided by qualified personnel on both a contractual or as-needed basis.

Varian, Inc. Analytical Instruments Sales Offices

For Sales or Service assistance and to order Parts and Supplies, contact your local Varian office.

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www.varianinc.com

Table of Contents

Introduction.....	1
Typographical conventions	1
Injection Principle	2
Full loop injections	2
Partial loopfill injections.....	5
µL Pick-up injections.....	6
Front view of the ProStar 430	10
Back view of the ProStar 430.....	11
Keypad and Display of the ProStar 430.....	12
Installation.....	15
Preparation for Use	16
HPLC connections	17
Waste tubing	17
Wash solvent.....	18
Syringe	19
Needle assembly.....	21
Optional Needles	23
Combination of Syringe, Sample Loop and Buffer Tubing	24
Plates and sample handling.....	27
Configuring the ProStar 430 for User Defined Plate Heights	28
Placing Reagent Vials/Transport Vials	30
Plate feeder.....	31
I/O Connections	35
TTL outputs	37
TTL inputs	39
Operation.....	43
Getting Started	43
Menus of the ProStar 430	43
Convenient working order	44
Types of Methods and links to Series.....	45
Executing a Series	46
Executing a Series in remote control	47
Using the ProStar 430.....	47
Ready Menu.....	55
System Menu	57

Methods Menu	62
Series Menu	70
Explanation	72
Programming Options	74
Programming Options	75
Maintenance and Troubleshooting	77
Analytical problems	78
Error Codes	82
Injection valve and ISS unit	82
Syringe dispenser unit	82
Injection needle unit	82
Plate	83
Vials	83
Electronics	83
Plate holder	84
Appendix	85
Specifications	85
General	85
Analytical performance	86
Programming	87
Physical	87
Electrical	87
Communication	88
Options	88
Options and Accessories	89

Figures

Figure 1 Front view of the ProStar 430.....	10
Figure 2 Back of the ProStar 430	11
Figure 3 Keypad and display	12
Figure 4 Wash solvent bottle	18
Figure 5 Replacing the syringe.....	20
Figure 6 Needle Assembly	22
Figure 7 Feeder Option	34
Figure 8 Contact Closures outputs	37
Figure 9 TTL output	39
Figure 10 TTL Input	41
Figure 11 ISS Valve, interconnected ports.....	66
Figure 12 Injection sequence with 3 calibration wells between every 5 wells.....	73
Figure 13 Fluidics Diagram.....	78

Introduction

The ProStar 430 is an AutoSampler that offers a wide variety of possibilities for HPLC. It incorporates among other things the following features:

- suitable for use with various types of plates by various manufacturers.
- capable of executing full loop injections, partial loopfill injections and μL pickup.
- accommodates use of various types of syringes (volumes of 100 μL , 250 μL , 500 μL , 1000 μL , or 10 mL).
- fully compatible with other laboratory equipment.
- apparatus designed for routine analysis and method development.

In addition options like plate cooling and a second switching valve are available for the ProStar 430. These options enhance the uses for the ProStar 430.

Typographical Conventions

The following typographical conventions are used in this manual:

Initial Capitals	are used for names of the menus: Ready Menu, Series Menu etc.
<CAPITALS>	are used to refer to soft function keys
bold type	is used for names of keys in the keypad
[0100]	is used to indicate values that have to be entered with the numerical keypad.

Injection Principle


The ProStar 430 offers three different methods of injection for an analytical run:

- | | |
|------------------|--|
| Full loop | The sample loop is completely (quantitatively) filled with sample resulting in extremely good reproducibility. |
| Partial Loopfill | The sample loop is partially filled with sample; this means low sample loss and programmable injection volumes. |
| μ L Pick-up | After aspiration of sample, the sample is transported into the loop with transport liquid (mobile phase); this means no sample loss. |


The ProStar 430 uses a syringe to aspirate the sample from a well into the sample loop. To prevent contamination of the syringe the ProStar 430 is equipped with a buffer tubing between the syringe and the injection valve. Wash solvent is used to remove the sample from the buffer tubing and sample needle, and to rinse the buffer tubing and sample needle.

The switching sequence for a full loop injection is schematically shown in the following figures:


Full loop Injections


The initial situation: the injector is in the INJECT position. The sample needle has entered the well after the air needle has pre-punctured the septum. Headspace pressure is applied through the outer air needle to ensure that no air or vapor bubbles are formed during sample aspiration.


The syringe dispenser aspirates the programmed flush volume from the sample well to fill the sample line with sample and remove wash solvents.

The injection valve switches in LOAD position, placing a "sharp" sample front at the inlet of the sample loop.

For full loop injections the sample loop is quantitatively filled by transporting two or more times the loop volume through the loop, depending on the volume of the loop.

The injection valve switches in INJECT position. The sample loop is now part of the HPLC mobile phase flow path: sample is transported to the column. The analysis time starts.

If one injection is to be done from each well or if a wash routine has to be performed after every injection, the needle withdraws from the well directly after the injection and, if programmed, directly performs a wash. After the analysis time a new sequence is started.


If more than one injection is done from the same well without a wash routine, the ProStar 430 withdraws a flush volume after the analysis time to compensate for diffusion of mobile phase from the rotor groove into the first part of the sample line during the analysis time. The flush volume between injections is not programmable and is always 50% of the programmed flush volume. If the total amount of sample withdrawn with the next injection from the well will exceed the total volume of the buffer tubing, the buffer tubing is emptied into the wash position before the next injection. The next fill sequence will then start with a full flush volume.

Air Segment

An air segment can be used to reduce the amount of flush volume. This air segment is at the front of the flush volume and will not be injected and therefore will not influence the injection. Use of an air segment can be enabled in the System Menu (General Menu).

With a standard needle the flush volumes must be: minimal 30 μL for injections with air segment and 35 μL for injections without air segment. If the samples are highly viscous it may be


necessary to program larger flush volumes and reduce the syringe speed for better performance.


A = with air segment; B = without air segment

Partial loopfill Injections

The switching sequence for a partial loopfill injection is schematically shown in the following figures:

For partial loopfill injections the sample loop is filled by transporting the programmed injection volume into the sample loop.


The injection valve switches into the INJECT position. The sample loop is now part of the HPLC mobile phase flow path: the sample is transported to the column. The analysis time starts.

The next injection sequence will start with a flush of 50% of the programmed flush volume, in case an injection from the same vial and no wash routine is programmed. Otherwise it will start with a flush of the programmed flush volume. If the aspiration of sample for the next injection will exceed the total volume of the sample buffer tubing, the buffer tubing is emptied before the next injection. The next injection will start with the programmed flush, see also the full loop injections.

Air Segment


An air segment can be used to reduce the amount of flush volume. This air segment is at the front of the flush volume and will not be injected. Use of an air segment can be enabled in the System Menu (General Menu).


A = with air segment; B = without air segment

μ L Pick-up Injections


The switching sequence for a μ L pick-up injection is schematically shown in the following figures.


The initial situation: the injection valve is in INJECT position. The sample needle has entered the vial of transport liquid (mobile phase, to avoid disturbance of the chromatogram with an additional peak of the transport solvent) after the air needle has pre-punctured the septum. The headspace pressure, applied through the outer air needle, ensures that no air or vapor bubbles are formed during wash solvent aspiration.


For the first injection after a wash or after emptying of the buffer tubing, the syringe dispenser aspirates transport liquid from the transport vial to fill the sample line with transport liquid and remove wash solvent.


The needle moves from the transport vial to the sample well. The injection valve is switched to the LOAD position.

The programmed injection volume is aspirated from the sample well.

The sample needle moves back to the transport vial. The sample is quantitatively transported into the loop, with transport liquid (mobile phase) from the transport vial.


The injection valve is switched to INJECT. The sample loop is now part of the HPLC mobile phase flow path: sample is transported to the column. The analysis time starts to run.

The next sequence will skip the first withdrawal of transport solvent, unless a wash routine is performed or the ProStar 430 has emptied the buffer tubing to waste. In those cases the sequence is completely repeated.

Air Segment

If an air segment has been programmed, it appears at the front of the first plug of transport liquid and at the front of every sample plug. Use of an air segment can be enabled in the System Menu (General Menu).

NOTE: The air segment at the front of the sample plug is injected into the HPLC system.

μL pickup injections with (A) or without (B) air segment.

NOTE: In case of μL pick-up injections there will be no air pressure (headspace pressure) on the sample well to prevent errors due to air expansion during switching from sample well to transport vial.

You are advised to switch off use of an Air segment in the System Menu if the μL pick-up injection method is used.

Front View of the ProStar 430

1. Cover (in open position)
2. Buffer tubing
3. Injection valve (Rheodyne)
4. Tubing holder
5. Needle unit
6. Screw to fasten cover
7. Transport/reagent vials
8. Plate
9. Plate holder
10. Drain wash position
11. Condensed water and leakage
12. Wash position
13. Wash solvent bottle
14. Syringe
15. Syringe waste tubing


Figure 1 Front View of the ProStar 430

Back View of the ProStar 430

1. Communication interface connector (RS232 or RS422/485)
2. CE-mark
3. I/O connectors
4. Mains switch
5. Mains input
6. Fuses and voltage selector
7. Type label
8. Fan (if plate cooling is installed)
9. Ventilation holes


Figure 2 Back of the ProStar 430

Keypad and Display of the ProStar 430


Figure 3 Keypad and display

Soft Function Keys

The label assigned to these keys depends on the menu that is active. The label of each key is shown in the bottom line of the display.

Cursor Keys

Used to move to a different field in the display, to move to a different field in a menu, or to make a displayed value higher or lower.

Numeric Keys

- **0 to 9:** to enter numerals in the various programming fields.
- **CL:** to clear a value in a field or replace it by NONE or AUTO.
- **E:** stands for Enter; to go through menu lines or to confirm a choice made in a menu or a value entered. The entered value is checked for validity and then saved.

Function Keys

Run Control Keys

Start/Stop: to start or stop automatic processing, or to reset the system after an error has occurred.

Hold/Cont.: to hold or continue the analysis time. The analysis time is extended by the period that Hold is active.

Interrupt: not used

Priority: to stop a run to process a priority sample before analyzing the rest of the programmed sample series. Before the run is interrupted processing of the present sample will be finished. As soon as the priority sample has been analyzed, the analytical run is resumed. A priority sample is a series of one well with an injection method, a wash method and a time base method defined in a template (only possible if the correct settings are entered in the System Menu).

Programming Keys

Series: to enter the Series Menu in which series can be defined for an analytical run.

Methods: to enter the Methods Menu in which methods can be programmed for use in an analytical run.

Menu: this key can only be used if [MENU] or [MN] is shown in the top right hand corner of the display. If pressed more possibilities of the menu are displayed.

System: to enter the System Menu in which system settings can be entered

General Keys

Escape: allows the user to leave the programming mode or go to a previous level in the menu. Entered values are checked for validity and then saved.

Help: to display help information; available only for a limited number of functions!

Installation

CAUTION

Do not install the ProStar 430 in areas subject to shock, dust, or in direct sunlight. Do not place it near a source of heat, as this will disturb tray cooling (if this option is installed).

The ProStar 430 must only be connected to power sources and apparatus with protective grounding.

Execute the following steps to install the ProStar 430:

1. Lift the ProStar 430 from the shipping container. Make sure you keep the apparatus upright; lift it by placing your hands under the ProStar 430.
2. Place the ProStar 430 in its operating location, preferably on the left hand side of the HPLC system. Make sure the ventilation holes are not obstructed. Allow the instrument to acclimatize for 1 hour.
3. Install plate holder, Figure 1 (5), in the ProStar 430 (on the flat surface underneath wash solvent bottle; as far to the left and to the back as possible).
4. Check whether local voltage matches voltage indicated on back panel of the ProStar 430.
5. Connect the power cord to the ProStar 430, Figure 2,(5).
6. Switch the ProStar 430 on by using the switch at the back of the apparatus, Figure 2 (4).

The ProStar 430 now starts up. The display indicates that a self-test and initialization procedure are executed. After completion of this procedure the following appears on the display:

10:26 TUESDAY DECEMBER 16-97 [MENU] READY (x.xx)
<EXCHANGE> <WASH> <SYR END> <UTILS>

Preparation for Use

This section describes procedures for replacement or installation of parts, either for maintenance or for preparation of the ProStar 430 for an analytical run.

NOTE: Please note that the power should **not** be switched off for any of the procedures described below as soft function keys in the Ready Menu must be used.

A number of items required for use of the ProStar 430 are factory-installed. However, after the installation procedure has been executed, check that the following have been correctly installed before you start to use the ProStar 430:

- HPLC pump & column connections
- waste tubing
- wash solvent bottle
- syringe
- needle assembly
- syringe, sample loop and buffer tubing
- plate holder and plate
- reagent vials/transport vials.

NOTE: Loosen the screw at the right-hand side of the cover, Figure 1 (6), and lift the cover to be able to execute the procedures described in this section.

HPLC Connections

Make the following connections:

- HPLC pump to port 1 of the injection valve
- HPLC column to port 6 of the injection valve

The instrument has been flushed with isopropanol before dispatch from the factory. Make sure that the mobile phase of your HPLC system is miscible with isopropanol, or start up with an intermediate solvent as mobile phase (disconnect the HPLC column).

NOTE:	It is very important that the contents of the sample loop are injected in back flush onto the column, therefore: do not exchange column and pump connections at the injection valve.
-------	--

Waste Tubing

Make the following connections:

Syringe waste: put the end of the syringe waste tube, Figure 1, (15) in a bottle placed next to the ProStar 430.

Drain wash-position: connect the hose to the drain wash connector of the ProStar 430, see Figure 1 (10); place the other end of the hose in a bottle placed on the floor. Through this drain all liquid dispensed to waste at the back of the plate is removed.

Condensed water and leakage: connect the hose to the drain port of the ProStar 430, see Figure 1 (11). Place the other end of the hose in a waste container on the floor. Through this hose all leakage solvents and condensed water (if Peltier cooling is installed) are drained.

Make sure that the flow path of the hoses is not obstructed in any way.

Wash Solvent

The ProStar 430 has a 250 mL wash solvent bottle. Execute the following steps to install the 250 mL wash solvent bottle:


Figure 4 Wash Solvent Bottle

1. Fill the wash solvent bottle with the appropriate wash solvent. Use of methanol (organic solvent, no buffers) or a mixture of water and isopropanol (80/20) is recommended. Before using the wash solvent, degas it with Helium or an ultrasonic bath.
2. Screw the bottle to the cap in the holder.
3. Place the holder in the ProStar 430 as indicated in Figure 4.
4. Put the wash solvent tube in the wash solvent.
5. Lower the cover of the ProStar 430
6. Select soft function key **<SYR END>** in the Ready Menu to fill the syringe.

7. Place the syringe back in home position by selecting soft function key **<SYR HOME>**.
8. Repeat steps 6 and 7 until the wash solvent tube and the syringe are completely filled.
9. Select soft function key **<WASH>** to perform a standard wash routine.
10. If any air remains in the syringe, select **<SYR END>** again to fill the syringe with wash solvent; press **<SYR HOME>** again to move contents to waste. Repeat if there is still air in the syringe and gently tap the syringe as wash solvent is dispensed to waste.

If you use an application that requires more than 250 mL of wash solvent for a complete run, install a longer tube (with flanged end for valve fitting) and place a larger bottle next to the ProStar 430. To fill the wash solvent tube, you may have to repeat the above-mentioned filling procedure (steps 6 and 7) a few times.

Syringe

The ProStar 430 is supplied with a 250 μ L syringe. It is also possible to use the ProStar 430 with a 100 μ L, 500 μ L, 1000 μ L, or 10 mL syringe. Execute the following steps to install a syringe.


Figure 5 Replacing the Syringe

1. Select soft function key **<SYR END>** in the Ready Menu to move the syringe to end position.
2. Lift the cover.
3. Unscrew the top of the syringe (turn clockwise).
4. Pull the bottom of the syringe towards you; you can now remove the syringe (refer to figure above).
5. Fill the new syringe with wash solvent and make sure that all air bubbles are removed from the syringe.
6. Connect the bottom of the filled syringe to the ProStar 430.
7. Screw the top of the filled syringe to the ProStar 430 (counterclockwise).

8. Lower the cover.
9. Select soft function key **<SYR HOME>** to remove air from the syringe. The syringe moves to home position and its contents is dispensed to waste.
10. If any air remains in the syringe, select **<SYR END>** again to fill the syringe with wash solvent; press **<SYR HOME>** again to move contents to waste. Repeat if there is still air in the syringe and gently tap the syringe as wash solvent is dispensed to waste.
11. Select soft function key **<WASH>** in the Ready Menu to execute a standard wash routine. All tubing connected to the syringe valve is filled and rinsed.


Select the appropriate syringe in the System Menu and install the correct buffer tubing.

Needle Assembly

The needle used for sampling consists of two parts:

- pre-puncturing needle: a hollow needle used for puncturing of the septum, capmat or sealer; also used to put headspace pressure on the sample (approximately 0.5 bar).
- sample needle: placed inside the hollow prepuncturing needle; used for the actual transport of sample. Different types of needles can be used here (see Appendix). If a needle with deviating diameter is used, a different air outlet nut, see Figure 6 (6), must be used that matches the injection needle.

NOTE: Most commercially available sealers or capmats cannot be used in combination with headspace pressure. You are advised to switch off headspace pressure in those cases (System Menu, General Menu).

The following elements can be identified in the needle assembly:

1. Nut and ferrule
2. Needle tubing
3. Needle connection nut
4. Standard sample needle
5. Needle holder
6. Air outlet nut
7. Prepuncturing needle
8. Sensor (detects plate height, missing plate and missing vials).

Numbers 1, 2, 3, 4, and 6 constitute the sample needle.

Figure 6 Needle Assembly

Execute the following steps to replace a needle:

1. Loosen the needle connection nut (3).
2. Loosen the nut (1) .
3. Carefully pull out sample needle and tubing.
4. Insert a new sample needle and tube through the needle holder (5) and tighten the nut (4).
5. Connect the other end of the tube to port 4 of the injection valve using a nut and ferrule. Do not tighten too much to prevent block of tubing.

6. Lower the cover of the ProStar 430.
7. Check sample needle height (default height: 5 mm from plate). If necessary, adjust the value in the System Menu (General Menu).
8. Select soft function key **<WASH>** in the Ready Menu to clean the new sample needle.

Optional Needles

Needles of the following types can be installed in the ProStar 430:

- LSV needle sample needle with large inner diameter for viscous samples, or in case large volumes are loaded in the loop
- PEEK needle biocompatible needle
- Fused silica needle needle with small inner diameter for small injection volumes
- Extended needle needle for switching valve placed in the side panel of the ProStar 430
- Silica-coated needle needle for use with liquids that would be damaging to a stainless steel needle

Refer to the Appendix for an overview of options available for the ProStar 430.

Execute the following steps to install any of these needles:

1. Remove needle as described in *Needle Assembly Section*, page 21.
2. Remove the standard air outlet nut, Figure 6 (6) and replace it by the nut supplied with the optional sample needle.
3. Install optional sample needle as described in *Needle Assembly* section, page 21
4. Adjust settings in System Menu (General Menu) to the volume of the new needle tubing, page 57).

Combination of Syringe, Sample Loop and Buffer Tubing

The 250 µL syringe is the standard syringe; combined with the standard 500 µL buffer and the standard 100 µL sample loop, the following injection volume range is available for the various injection modes:

Full loop	100 µL
Partial loopfill	1 - 50 µL
µL pick-up	1 - 27 µL

The maximum injection volumes are calculated with the following formulas:

Full loop:	injection volume	= loop volume
Partial loopfill:	max. injection volume	= 50% of loop volume
µL pick-up:	max. injection volume	= (loop volume - 3 x needle tubing volume) 1/2

Five sizes of syringes can be used in the ProStar 430: 100 µL, 250 µL, 500 µL, 1000 µL and 10 mL. The 10 mL syringe can only be used in the User Program (must be enabled in System Menu) and not for the standard injection modes.

Full loop gives maximum reproducibility (RSD < 0.3%), but not maximum accuracy, since loop volume is specified with an accuracy of ± 10%.

Minimum sample loss = 230 µL (2 x loop overfill + flush volume for needle) for the standard 100 µL loop.

Partial loopfill gives maximum accuracy (depends on syringe accuracy) and reproducibility better than 0.5% RSD for injection volumes > 5 µL

Minimum sample loss (Flush volume) = 30 µL

30 μ L is the recommended minimum flush volume (combined with an air segment); smaller flush volumes can be programmed, but will result in decreased performance.

μ L Pick-up means zero sample loss, maximum accuracy (same as partial loopfill), but slightly diminished reproducibility: RSD better than 1% for injection volumes > 5 μ L.

5 μ L of air is injected together with the sample, if an air segment is selected in the System Menu.

For some cases other combinations of syringe, loop and/or buffer are advised:

Injection volumes smaller than 5 μ L:

Partial loopfill: use a 100 μ L syringe for maximum reproducibility and accuracy. Use a 20 μ L sample loop to avoid loss of accuracy due to expansion of the loop content when switching from inject to load position prior to sample loading. When working with high pressure (200 bar), this loss may be 0.1 - 0.5 μ L for a 100 μ L loop.

Note that the minimum sample loss in partial loopfill mode is 30 μ L (recommended minimum flush volume) for the first injection and an additional 15 (always half the programmed flush volume) for additional injections from the same well. If a wash between injections has been programmed, sample loss is 30 μ L for every injection. For zero sample loss injections, use the μ L Pick-up injection mode.

μ L Pick-up: use a 100 μ L syringe for optimum accuracy and reproducibility. Do not use a smaller sample loop! The sample plug is transported into the loop, preceded by a programmable air segment of 5 μ L (if selected in System Menu; you are advised to switch off the air segment with μ L Pick-up), with a plug of transport liquid which equals 2.5 times the programmed needle tubing volume.

Injection volumes up to twice the standard:

With the standard 250 μL syringe, standard needle with tubing (15 μL) and standard 500 μL buffer, but with a 200 μL sample loop, the maximum injection volumes are:

Full loop	200 μL (sample loss remains 230 μL since loops > 100 μL need only one loop volume overfill; 30 μL pre-flush)
Partial loopfill	100 μL
μL Pick-up	77 μL

Volumes smaller than 5 μL may be injected, but reproducibility and accuracy may not be < 0.5% for partial loopfill or < 1% for μL pick-up.

In short: loop < 100 μL : loop has to be filled three times

loop \geq 100 μL - 499 μL : loop has to be filled twice

loop \geq 500 μL : loop has to be filled 1.5 times

For volumes larger than 200 μL :

Use the 2000 μL buffer, the appropriate sample loop size and the appropriate syringe: Syringe volume > 2 x injection volume. Injection volumes larger than 500 μL are possible, but the sample may contaminate the syringe. Program sufficient wash after use!

Appropriate buffer tubing for each type of syringe:

Syringe	Buffer tubing
100 μL	500 μL
250 μL	500 μL
500 μL	2000 μL
1000 μL	2000 μL
10 mL	15 mL

NOTE: The 10 mL syringe and the 15 mL buffer tubing can only be used in combination with the User Program and not with the standard injection methods.

Plates and Sample Handling

The ProStar 430 accommodates use of the following types of plates:

- 96-low wells
- 96-high wells
- 384-low wells
- 48-vials
- Plates of user defined height (These must be specially configured.)

A sensor, Figure 6 (8), takes care of plate detection, plate height detection, and vial detection.

Because the ProStar 430 uses headspace pressure during sample injections, it is very important that samples are properly handled. Note the following:

- Standard wells can best be filled by means of a narrow-end pipette to allow air to escape when filling the well.
- Do not fill wells to the rim as sample might be forced into the pre-puncturing needle, causing cross-contamination of samples and contamination of the sample needle.
- For vials: make sure the seals are airtight to prevent air bubbles in the sample and prevent evaporation of volatile samples; check seals after crimping: if the cap can be turned easily, the seal is not airtight and the handcrimper should be adjusted.

NOTE: If wells are used that are not airtight switch off headspace pressure in the System Menu (General Menu).
Check whether the needle height is sufficient for the new type of plate that is installed (System Menu, General Menu). Adjust if necessary.

To replace a plate in the ProStar 430:

1. Select soft function key **<PLATES>**.
2. Select soft function key **<EXCHANGE>** in the Ready Menu. The plate moves to the left.
3. Take out and replace the plate.
4. Select soft function key **<PLATE HOME>**. The plate moves to operating position again.

If you have replaced the plate by a plate of the same type, you are now ready. If you have installed a new type of plate, execute the following steps:

1. Press **System**.
2. Select soft function key **<PLATES>**.
3. Press **E**.
4. Select the soft function key for the type of plate concerned.
5. Press **E** and determine whether to process the plate in **<ROWS>** or **<COLUMNS>**.
6. Press **E**.
7. Press **Escape** twice to return to the Ready Menu.

A message appears to indicate that all programmed series will be reset. The user must redefine series because the settings in the System Menu have been altered.

Configuring the ProStar 430 for User Defined Plate Heights

It is possible to use plates of almost any height in the ProStar 430. In order to do this, the ProStar 430 must have version 2.02 or later software. The system must be configured to use user

defined plates, and it must be calibrated for the actual plate height.

Execute the following steps for configuration:

1. Press **System**
2. Press **<PLATES>**.
3. Press **E**.
4. Press **Menu** The menu entry user defined will be displayed.
5. Press **<User Def>**
6. Press **E**.
7. Enter the number of plates that will be used in the instrument. Note, the Maximum number of plates is equal to the number of plates in the plate holder + 1 (the one that is being sampled at the present time.)
8. Press **E**.
9. Enter whether the plates will be processed in rows or columns and whether transport vials will be used.

The ProStar 430 is now configured to have a user defined plate height. Now it will be necessary to calibrate the plate height.

To calibrate the plate height:

1. Press **Escape** until the Ready menu is accessed
2. Press **<Plates>** to Access the plates menu
3. Press **<Plate>** to set the plate height

A screen will appear on which you can press **Help** or **E**. If you press Help, the following text will appear on multiple screens to instruct you how to adjust the plate height:

*Press **<Default>** to accept the following heights
As well for single as for the maximum number
Of stacked plates.*

*One – 96 – low plate 14.5 mm
One – 384 – Low plate 14.5 mm
One – 96 – high plate 42.0 mm
Six stacked – 96 – low plates 77 mm*

*Six stacked - 384 - low plates 79 mm
Two Stacked – 96 – high plates 82 mm*

*Press **<ADJUST>** to calibrate different plate
Heights*

*Convention in all adjust screens:
Press **E** to accept new plate height
Press **ESC** to exit without saving*

*Place maximum number of plates in lift
This number depends on plate settings!
Place a ruler on top the right side
Press **<LIFT-DOWN>** or **<LIFT-UP>** until the top
Plate touches the ruler
Press **E** to enter the next adjust screen*

*Lift will move to maximum plate height
Take out plates. Leave 1 plate in the lift
Place a ruler on top of right lift again*

*Press **<LIFT-DOWN>** or **<LIFT-UP>** until
The top plate touches the ruler
Press **E** to accept single plate height*

*Lift will move to home position
Screen will return to main menu*

Placing Reagent Vials/Transport Vials

Execute the following steps to replace reagent vials/transport vials:

1. Select soft function key **<EXCHANGE>** in the Ready Menu. The plate holder moves to the left.
2. Take out the reagent vials/transport vials, Figure 1 (7), and replace them with other reagent vials/transport vials.

3. Select soft function key **<PLATE HOME>**. The plate moves to operating position again.

NOTE: Reagent and transport vials can be placed in any of the four positions. Transport vials must be placed in a continuous row.

In the System Menu, go to the Plates Menu to define the position for first and last Transport vial.

Plate Feeder

The Plate feeder is a factory-installed option that allows you to feed a maximum of 7 plates of the same type into the Endurance (one plate on the plate holder and six in the feeder). The following types of plates can be used with the feeder:

- 96-low plates (max. 6 plates in feeder)
- 96-high plates (max. 2 plates in feeder)
- 384-low plates (max. 6 plates in feeder)

The maximum height of plates that can be stacked into the feeder lift is 83.5 mm.

CAUTION

The 48-vials cannot be used in the plate feeder.

Changes in operation if the Endurance Feeder option is installed:

In the Series menu you have to enter a plate number to define the location of wells, apart from the location of wells in rows and columns. For example, you can define the location of sample wells in a mix method on plate 1, and define the location of the destination wells on plate 5.

The Ready menu displays <PLATES> as the first soft function key; press this key to enter the Direct control plate handling menu. This menu offers the following possibilities:

- <EXCHANGE>** to move the plate holder (Figure 1, number 9)
- <LEFT-LIFT>** to move the left lift up or down (refer to Figure 7)
- <RIGHT-LIFT>** to move the right lift up or down (refer to Figure 7)
- <PLATE>** To enter the plate height determination menu. This menu offers the following possibilities

<DEFAULT>

For default plate sizes:

Height of one 96-low plate: 14.5 mm

Height of six stacked 96-low plates: 77 mm

Height of one 96-high plate: 42 mm

Height of two stacked 96-high plates: 82 mm

Height of one 384-low plate: 14.5 mm

Height of six stacked 384-low plates: 79 mm

<ADJUST>

Before adjusting the plate height, be sure the appropriate plate type has been selected in system settings (System Menu, Plates, refer to page 47).

The Adjust Menu allows you to calibrate the Endurance feeder to use plates of a different height:

1. Place six or two plates (depending on settings) in the right lift (refer to Figure 7).
2. Place a ruler on top of the right lift
3. Select **<LIFT-DOWN>** or **<LIFT-UP>** until the top plate touches the ruler.
4. Press **E** to confirm the height of the six stacked plates.
5. Remove 5 plates from the right lift.
6. Repeat steps 2 – 5 to calibrate the Endurance for single plate height.
7. Press **E** to confirm the new plate height settings.

Using the Feeder:

1. Make sure there are no plates in the left lift.
2. Define the type of plates you are going to use in the usual manner (refer to page 27).
3. Select **<PLATES>** in the Ready menu.
4. Select **<RIGHT-LIFT>** in the Plates menu.
5. Select **<LIFT-UP>**. The lift moves up.
6. Place a plate on the right lift (see Figure 7, number 2).
7. Select **<NEXT PLATE>** to move the lift one position down. You can now place another plate on the right lift.
8. Repeat step number 7 until:
 - you have placed the required number of plates on the lift; press **Escape** twice to return to the Ready Menu.
 - until the lift is full; press **Escape** to return to the Ready Menu.

Keep in mind that the plate that is on the plate holder in the Endurance is number 1, the plates on the lift are numbered 2 through 7. The uppermost plate in the lift is number 2.

The Feeder is now ready for use.

The left lift is used for storage of plates that have been processed and for plate switching sequences. At the end of the programmed series Endurance will make sure that plate number 1 is always on bottom of the left lift, and the plate with the highest number will be on the plate holder.

Example: If 7 plates are used, plate number 7 will be on the plate holder at the end of the programmed series, plate number 1 will be at the bottom of the left lift, plate number 2 will be on top of number 1 and the uppermost plate in the left lift will be number 6.

CAUTION

Make sure that the tray holder of the Endurance always contains a plate at the start of programming a series.


Figure 7 Feeder Option

- | | |
|-------------------------|--------------------------------|
| 1. Grab head | 4. Drain wash-position |
| 2. Right lift | 5. Condensed water and leakage |
| 3. Plates in right lift | 6. Left lift |

I/O Connections

The ProStar 430 has six standard I/O connectors on the rear side; five OUTPUT connectors and one INPUT connector. Refer to Figure 2 for location of connectors.

The communication connector is a standard RS232 or an RS422/485 communication interface connector. The configuration of the I/O connector is described in the tables below.

CAUTION

The manufacturer will not accept any liability for damages directly or indirectly caused by connecting the ProStar 430 to instruments which do not meet relevant safety standards.

Contact closures outputs

The following three tables with programmable outputs (P1), marker outputs (P4) and auxiliary outputs (P5) are contact closures outputs (floating NO/NC contact); see Figure 8.

Table 1 Connector P1 OUTPUTS (2 programmable outputs and alarm output)

1. OUT 1 - Normally open	8. Spare
2. OUT 1 - Common	9. Spare
3. OUT 1 - Normally closed	10. Alarm output - Normally open
4. OUT 2 - Normally open	11. Alarm output - Common
5. OUT 2 - Common	12. 24 V DC
6. OUT 2 - Normally closed	13. Power ground
7. Spare	14. Power ground

$$V_{MAX} = 28 V_{DC} / V_{AC}, I_{MAX} = 0.25 A$$

NOTE: The Alarm output will be activated whenever an error occurs. Refer to the *Maintenance and Troubleshooting* section on page 82 for a description of the error codes of the ProStar 430.

Table 2 Connector P4 MARKERS

1. Inject marker - Normally open	8. Labeled well marker - Common
2. Inject marker - Common	9. Labeled well marker - Normally closed
3. Inject marker - Normally closed	10. STOP I/O - Normally open
4. Well marker - Normally open	11. STOP I/O - Common
5. Well marker - Common	12. STOP I/O - Normally closed
6. Well marker - Normally closed	13. 24 V DC
7. Labeled well marker - Normally open	14. Power ground
	15. Power ground

$$V_{MAX} = 28 V_{DC} / V_{AC}, I_{MAX} = 0.25 A$$

Table 3 Connector P5 AUXILIARIES

- | | |
|----------------------------|-----------------------------|
| 1. AUX 1 - Normally open | 8. AUX 3 - Common |
| 2. AUX 1 - Common | 9. AUX 3 - Normally closed |
| 3. AUX 1 - Normally closed | 10. AUX 4 - Normally open |
| 4. AUX 2 - Normally open | 11. AUX 4 - Common |
| 5. AUX 2 - Common | 12. AUX 4 - Normally closed |
| 6. AUX 2 - Normally closed | 13. 24 Vdc |
| 7. AUX 3 - Normally open | 14. Power ground |
| | 15. Power ground |

$V_{MAX} = 28 \text{ Vdc} / \text{Vac}$ $I_{MAX} = 0.25\text{A}$
--

NOTE: Maximum current for 24 Vdc supply is 0.5A total.


Figure 8 Contact Closures outputs

TTL Outputs

The following tables show the marker outputs (P2) and a 4 bit time base code output (P3), programmable in a time base method. Both connectors are TTL level outputs, a connection diagram is shown in Figure 9.

Table 4 Connector P2 TTL OUTPUTS

1. INJECT MARKER	9. not connected
2. VIAL/WELL MARKER	10. not connected
3. LABELED WELL MARKER	11. not connected
4. STOP I/O	12. not connected
5. not connected	13. Signal ground
6. not connected	14. Signal ground
7. not connected	15. Signal ground
8. not connected	

All markers are active low (logical 0).

$V_{MAX} = 5.5 \text{ V}$, logical 1 > 3.5 V, logical 0 < 1.0 V. DC output source / sink current $\pm 20 \text{ mA}$.

NOTE: A marker output pulse will be generated when the injection valve switches from LOAD to INJECT. However, in a User Program markers have to be programmed by the user.

Table 5 Connector P3 TIMED OUTPUTS; 4 bit time base code output

1	TB 0 (HEX) (1) 6	Signal ground
2	TB 1 (HEX) (2) 7	Signal ground
3	TB 2 (HEX) (4) 8	Signal ground
4	TB 3 (HEX) (8) 9	Signal ground
5	not used	

$V_{MAX} = 5.5 \text{ V}$, logical 1 > 3.5V, logical 0 < 1.0V. DC output source / sink current $\pm 20 \text{ mA}$.


Figure 9 TTL output

TTL Inputs

The following connector is an active high or active low TTL input; it can be defined in the System Menu. The NEXT INJECTION INPUT and the NEXT WELL INPUT can be used when the ProStar 430 works in REMOTE CONTROL. The FREEZE INPUT and STOP I/O input can be used to control the ProStar 430 by other devices. The four inputs (INPUT 1 to 4) can only be used in the user program, e.g. to control the sequence of the steps in this method. A connection diagram is shown in Figure 9.

Table 6 Connector P6 INPUTS (TTL)

1. NEXT INJECTION INPUT	9. Signal ground
2. NEXT WELL INPUT	10. Signal ground
3. FREEZE INPUT	11. Signal ground
4. STOP I/O	12. Signal ground
5. INPUT 1	13. Signal ground
6. INPUT 2	14. Signal ground
7. INPUT 3	15. Signal ground
8. INPUT 4	

Next Injection Input This input will start the next injection sequence when the ProStar 430 is started in remote control. When the injection sequence is finished the ProStar 430 will wait for the next input.

From the Ready Menu a NEXT INJECTION INPUT will start the last programmed series. In this case the ProStar 430 will not wait for the NEXT INJECTION INPUT before continuing with the next injection. The ProStar 430 will execute the complete RUN as if it was started with the **Start/Stop** key.

Next Well Input With this input the ProStar 430 will perform the next injection from the next well, even if not all injections from that well in the programmed injection method have been executed.

Freeze Input The ProStar 430 will freeze the analysis time for the time this input is active. If the FREEZE INPUT is activated while the analysis time is not running, the ProStar 430 will perform all programmed pre-injection sample handling (mix method and loading part of the injection method). But the ProStar 430 will wait with injecting the sample until the FREEZE INPUT is no longer active.

Stop I/O With this input the run of the ProStar 430 is immediately aborted. The Ready Menu appears in the display. In case the ProStar 430 is in remote control, the run of the ProStar 430 is immediately aborted but the ProStar 430 remains in remote control and cannot be restarted with a NEXT INJECTION INPUT

Input 1-4

Programmable input, can be used in the user program.


Figure 10 TTL Input

Operation

Getting Started

Menus of the ProStar 430

The software of the ProStar 430 is menu-driven. The most important menus are:

- **Ready Menu:** this menu appears after the ProStar 430 has been switched on. It offers very general possibilities for among other things methods management and communication with other apparatus.
- **System Menu:** this menu appears after the **System** key in the keypad has been selected. It offers general possibilities for entering settings for an analytical run. Please note that settings may have been factory-installed to suit particular wishes of the user. In that case no changes will have to be made in the System Menu. You are advised to enable only those facilities that you will need, to ensure that menus are as concise as possible.
- **Methods Menu:** this menu appears after the **Methods** key in the keypad has been selected. The menu allows you to program a method to be used in an analytical run and to assign a number to it.
- **Series Menu:** this menu appears after the **Series** key in the keypad has been selected. The menu allows you to define a series, and to assign a number and a method to it for the analytical run.

Use the following **keys** to go through the menus:

E	ENTER; to confirm a choice or to select a choice made in a screen, or to step through menu lines
cursor keys	to change values in a field or to move to a different field in a screen
Escape	to return to a previous menu
CL	to remove a value from a field and enter NONE or DEFAULT
soft function keys	labels displayed in bottom line of the screen; to go to submenus

If [MENU] or [MN] is displayed in the top right hand corner of the screen, you can press the Menu key in the keypad to display more possibilities offered by the menu. Refer to page 12 for an explanation of all keys in the keypad.

Convenient Working Order

After you have determined what type of analytical run you want to perform, the most convenient working order for the ProStar 430 is:

1. Enter settings in the System Menu. Note that the settings for the ProStar 430 probably already have been correctly entered (factory-installed).
2. Program a method for the analyses you wish to perform in the Methods Menu.
3. Define a series and link a programmed method to a range of wells in the Series Menu.
4. Execute the series.

Please note that you may use a different order; however, you should keep in mind that settings entered in the System Menu determine which possibilities appear in the Methods and Series menus.

Types of Methods and Links to Series

The ProStar 430 offers the following types of methods for different parts of the sample handling routine:

- **injection** method: contains information on the injection routine, flush volume and analysis time.
- **wash** method: describes a wash volume and when a wash must be executed.
- **mix** method: a pre-injection method in which additional sample handling can be performed (e.g. pre-column derivatization).
- **timebase** method: a post-injection method with which outputs to other devices (e.g. integrator or pump) and switching of the ISS valve are controlled.
- **user program**: offers the possibility to program sequences of all actions that can be executed by the ProStar 430 in separate steps.

Each programmed method is assigned a number. The ProStar 430 offers the possibility to store a combination of defined methods in a template. A template is also identified by a number.

Methods must be linked to series before they can be used. The following possibilities are offered by the ProStar 430:

- You can assign an individual method to a series: methods (mix, injection, wash, timebase) can be linked to wells in a series.
- You can assign a template to a series: a combination of various programmed methods (mix, injection, wash, timebase) can be defined in a template. The template is linked to a range of wells in a series. In this way all steps in an analytical run are laid down and stored.
- You can assign a user program to a series: a possibility to combine all possible steps in the analytical process in

one program. The user determines the order of the separate actions the ProStar 430 has to perform.

Executing a Series

Execution of a series is only possible if you have programmed a method and defined a series for the samples you wish to analyze. Series are not stored in battery backup and exist only for as long as the ProStar 430 is switched on.

Execute the following steps:

1. Start programming the series by pressing **Start/Stop**.
2. Enter the number of the first series to perform and the number of the last series to perform.
3. Select **<START>** to start the actual analytical run. The ProStar 430 starts to execute the series you have defined.

After the ProStar 430 has executed the run, the Ready Menu will appear again.

Refer to *Using the ProStar 430*, page 47 for a number of examples illustrating this working order.

It is possible to program series and methods during a run. Press Series or Methods; the possibilities offered in the menus are identical to those offered when the ProStar 430 is idle.

If a series or method is changed, the new values become active the next time the ProStar 430 starts a series. The series currently running are not affected by the changes.

Executing a Series in Remote Control

To execute a series from remote control, execute the following steps:

1. Press **Start/Stop**.
2. Enter the number of the first and the last series to be performed.
3. Select **<REMOTE>** to enter the remote control mode. The ProStar 430 will now operate as slave of another device and can be controlled with Next injection input and Next well input. To indicate that remote control is active, an “r” is displayed in the bottom left corner of the display during execution of the series. At the end of the series the message “Series completed via remote control” is displayed.
4. Press **Escape** to return to the Ready Menu.

Refer to the I/O Connections section beginning on page 35 for more information on remote control.

Using the ProStar 430

Following are examples of actions that can be performed with the ProStar 430. Please note that this chapter does not describe all types of actions that can be performed. Try to do these examples to learn to work with the ProStar 430.

These examples can be executed after the ProStar 430 has been installed and all items described in the previous section correctly set up.

Example 1. A 10 μ L partial loopfill injection

After the ProStar 430 has been switched on, an initialization procedure is executed. In this example it is assumed that a loop of 100 μ L, needle tubing of 15 μ L and a syringe of 250 μ L has been installed. Wait until the Ready Menu appears on the screen, then execute the following steps:

For system settings:

Press keys	Description
System	to enter the System Menu
<GENERAL> E	to enter the General Menu
[0100] E	to define the volume of the installed loop
[015] E	to define the volume of the needle tubing
<250> E	to define the volume of the syringe
<NORMAL> E	to set syringe speed to normal
[05] E	to set sample needle height to 5 mm
<YES> E	to enable use of air segment
<NO> E	to switch off headspace pressure
Escape	to return to the System Menu
<PLATES> E	to enter the Plates Menu
<96-LOW>	to define the type of plate to be used
<IN ROWS>	to define the processing order of wells
Escape Escape	to return to the Ready Menu

For this example all other settings used will be default.

To program a method:

Press keys	Description
Methods	to enter the Methods Menu
<INJECTION> [01] E	program injection method number 1
<PARTIAL> E	to select partial loopfill injection method
[100] E	to define an analysis time of 1 minute
[030] E	to define a flush volume of 30 µL
[1] E	to define the number of injections per well
[100] E	to set the injection volume at 10.0 µL
Escape Escape	to return to the Ready Menu

To define the series:

Press keys	Description
Series	to enter the Series Menu
[01] E	to define the Series number
[01] E	to define the injection method number
CL E	to enter <NONE> for wash method
<ROW A> [01] E	to define location of the first sample well
<ROW A> [01] E	to define location of the last sample well
Escape	to return to the Ready Menu

To run the series:

Place a sample in position A 1 of the plate.

Press keys	Description
Start/Stop	to start the ProStar 430
[01] E	to start at series number 1
[01] E	to stop after execution of series number 1
<START>	to start the analytical run

The ProStar 430 will now locate well A 01 and perform a 10 µL partial loopfill injection. The display of the ProStar 430 will indicate the status (Checking tray, Flushing, Loopfill, Running, Rinse buffer, Running). The display also indicates the number of the defined series (01), the method number (01) and the well on which the analysis is performed (A 01).

At the end of the defined analysis time the Ready Menu will be displayed again to indicate that the ProStar 430 is ready for the next analytical run.

Example 2. A 3 x 10 µL injection with µL pick-up, wash between injections

In this example a different injection method is used than in the previous one. For that reason several settings must be adapted in the System Menu.

For system settings:

Press keys	Description
System	to enter the System Menu
<GENERAL> E	to enter the General Menu
E until Air segment appears	to go to the Air segment field
<NO> E	to switch off air segment
Escape	to return to the System Menu
<PLATES> E	to enter the Plates Menu
E E	to go to transport vials field
[1] E	to define position of the first transport vial
[1]	to define position of the last transport vial
Escape Escape	to return to the Ready Menu

To program a method:

Press keys	Description
Methods	to enter the Methods Menu
<INJECTION> [02] E	to define method number 02
<PICK-UP> E	to select the injection mode for this method
[100] E	to define the analysis time
[3] E	to define the number of injections per well
[100] E	to define volume of 10.0 µL for 1st injection
[100] E	to define volume of 10.0 µL for 2nd injection
[100] E	to define volume of 10.0 µL for 3rd injection
Escape	to return to the Methods Menu
<WASH>	to enter the Wash Menu
[01] E	to define wash method number 01
<INJECTION> E	to select wash between injections
[300]	to define the wash volume
Escape Escape	to return to the Ready Menu

To define the series:

Press keys	Description
Series	to enter the Series Menu
[01] E	to define the series number
[02] E	to define the injection method for this series
[01] E	to define the wash method for this series
<ROW A> [01] E	to define the location of the first sample well
<ROW A> [01]	to define the location of the last sample well
Escape	to return to Ready Menu

To run the series:

Put a vial with transport solvent (mobile phase) in transport vial position 1 (left). Make sure the transport vial is correctly filled before starting a new series.

Press keys	Description
Start/Stop	to start the ProStar 430
[01] E	to start at series 01
[01] E	to stop after series 01
<START>	to start execution of the series.

At the end of the defined analysis time the Ready Menu will appear again to indicate that the ProStar 430 is ready for the following next run.

Example 3. A 1:10 dilution followed by a 10 µL partial loopfill injection

This example describes how to let the ProStar 430 transfer 360 µL from Reagent A to the destination vial, add 40 µL of sample, mix 3 times with 250 µL and subsequently inject 10 µL.

For system settings:

Press keys	Description
System	to enter the System Menu
<USAGE> E	to enter the Usage Menu
E E	to go to the Mix field
<ENABLED>	to enable use of mix methods
Escape	to return to the System Menu
<PLATES> E	to enter the Plates Menu
<96-HIGH>	to select 96 deepwell plate
<COLUMNS>	to select processing in columns (A1, B1, etc)
CL	to disable use of transport vials
Escape Escape	to return to the Ready Menu

As soon as a change has been entered in the System settings, the message “ALL SERIES DEFAULT” appears. The user will have to redefine series because the settings have been changed.

To program the injection method:

Press keys	Description
Methods	to enter the Methods Menu
<INJECTION> [03] E	to enter the Injection Menu
<PARTIAL> E	to select partial loopfill injection mode
[100] E	to define the analysis time
[50] E	to define the flush volume
[3] E	to define the number of injections per well
[100] E	to enter the injection volume for 1st injection
[100] E	to enter the injection volume for 2nd injection
[100] E	to enter the injection volume for 3rd injection
Escape	to return to the Methods Menu

To program the mix method:

Press keys	Description
<MIX>	to enter the Mix Menu
[1] E	to define Mix method number 1
<INSERT>	to define mix method step number 1
<ASPIRATE> [50] <AIR> E	to aspirate an air segment of 5.0 µL
<INSERT>	to define mix method step number 2
<ASPIRATE> [2250]	to aspirate 225.0 µL
Menu <REAG-A> E	to perform defined asp. from reagent vial A
<INSERT>	to define mix method step number 3
<DISPENSE> [1800]	to dispense 180.0 µL to
▶ ▶ [02] E	destination well
<INSERT>	to define mix method step number 4
<DISPENSE> [500]	to dispense 50.0 µL
<WASTE> ▶ [5] E	to dispense defined volume to waste
<INSERT>	to define mix method step number 5
<REPEAT> [1] ◀ [4] E	to repeat last four steps once
<INSERT>	to define mix method step number 6
<ASPIRATE> [50] <AIR> E	to aspirate an air segment of 5.0 µL
<INSERT>	to define mix method step number 7
<ASPIRATE> [600] <SAMPLE>	to aspirate 60.0 µL of
E	sample
<INSERT>	to define mix method step number 8
<DISPENSE> [400]	to dispense 40.0 µL
<DESTINATION> ▶ ▶ [2] E	to destination well
<INSERT>	to define mix method step number 9
<DISPENSE> [500] <WASTE>	to dispense 50.0 µL
▶ [5] E	to waste
<INSERT>	to define mix method step number 10
<ASPIRATE> [500] <AIR> E	to aspirate an air segment of 50.0 µL
<INSERT>	to define mix method step number 11
<ASPIRATE> [2000]	to aspirate 200.0 µL
<DESTINATION> ▶ [3] E	from the destination well
<INSERT>	to define mix method step number 12
<DISPENSE> [2000] ▶ [9] E	to dispense 200.0 µL to the destination well
<INSERT>	to define mix method step number 13
<REPEAT> [3]	to repeat the last 2 steps three times
Escape Escape	to return to the Ready Menu

To define the series:

Press keys	Description
Series	to enter the Series Menu
[01] E	to define series number 1
[01] E	to select Mix method number 1 for this series
[03] E	to select Injection method number 3
CL E	to select <NONE> for wash method
<ROW A> [1] E	to define location of first sample well
<ROW A> [1] E	to define location of last sample well
<ROW B> [1] E	to define location of first destination well
[1] E	to define position of Reagent A
Escape	to return to the Ready Menu

To run the series:

Put sample in position A 01; position B 01 is used as empty destination well. Place a filled reagent vial in position 1. Make sure the reagent vial is filled correctly before starting a new series.

Press keys	Description
Start/Stop	to start the ProStar 430
[01] E	to start at series number 1
[01] E	to stop after series number 1
<START>	to start processing of sample

The ProStar 430 will now start searching for the Reagent vial and transport 180 µL to the destination well twice, then 40 µL of sample will be added and after mixing 3 times a 10 µL injection will be performed.

Example 4. Defining a template and adding a protection code

This examples describes how to incorporate the injection method (02) and wash method (01) defined in example 2 in a template. A protection code will be added.

For system settings:

Press keys	Description
System	to enter the System Menu
<USAGE> E	to enter the Usage Menu
[123456] E	to enter a 6-digit code (memorize this code!)
E	to go to the mix methods field
<DISABLED> E	to disable use of mix methods
E E	to go to the template field
<ENABLED>	to enable templates
Escape Escape	to return to the Ready Menu

After use of templates has been enabled the message "ALL SERIES DEFAULT" appears. The user will have to redefine series because the settings have been changed.

To select the methods to be incorporated in the template:

Press keys	Description
Methods	to enter the Methods Menu
[123456] E	to enter the methods protection code
<TEMPLATE>	to enter the Template Menu
[01] E	to define the number for the template
[02] E	to define the injection method for this template
[01]	to define the wash method for this template
Escape Escape	to return to the Ready Menu

To define the series:

Press keys	Description
Series	to enter the Series Menu
[01] E	to define the Series number
[01] E	to define the Template method number
<ROW A> [01] E	to define the first sample well
<ROW B> [01] E	to define the last sample well
Escape	to return to the Ready Menu

To run the series:

Press keys	Description
Start/Stop	to start the ProStar 430
[01] E	to start analysis at series 01
[01] E	to stop after analysis of series 01
<START>	to start the analytical run

The ProStar 430 now performs the same actions as in Example 2, except that analysis is performed on two wells: A 01 and B 01.

NOTE: Select **<DEFAULT ALL>** in the Ready Menu (Utilities Menu) to erase all series and methods defined in these examples and to default all settings.

Ready Menu

NOTE: This section describes all possibilities offered by the ProStar 430 software, in the order in which they appear on the screen.

The Ready Menu contains the following soft function keys:

- <PLATES>** Use this key to move the plate to the left, move the lifts up and down and to adjust the plate height
- <EXCHANGE>** Use this key to move the plate to the left. In this position the plate can be replaced without damage to the equipment. Select soft function key **<Plate Home>** to move the plate to the operating position again
- <LEFT-LIFT>** Use this key to move the left lift up to the top to more easily take plates out of the lift.
- <RIGHT-LIFT>** Use this key to move the right lift up to the top to more easily take plates out of the lift.
- <PLATE>** Use this key to adjust the plate height for different size plates. The **<help>** key contains complete instructions for making the adjustment.
- <WASH>** Use this key to start a standard wash procedure. All tubing connected to the syringe valve will be filled and rinsed with wash solvent.
- <SYR END>** Use this key to move the syringe to end position if you wish to replace the syringe needle or to simplify filling of wash solvent tubing. A syringe volume of wash solvent is aspirated from the wash solvent bottle and the wash solvent tube is filled.

Select soft function key **<SYR HOME>** to dispense the syringe contents to syringe waste and to move the syringe to standard operating position again
- <UTILS>** Use this key to go to the Utilities Menu. If use of a method protection code is enabled in the System Menu, the code must be entered to access the Utilities Menu. The menu offers the following possibilities:

- <COPY>** to copy a method. Enter the type (mix, injection, timebase, wash) and the number of the method to be copied. Then enter a number to define the destination method. Any existing method stored under that number will be overwritten.
- <ERASE>** to erase a method (Template, Methods, User Program). If Template and User Program are disabled in the System Menu, the soft function keys for erasing a standard Method (mix, injection, wash, timebase) appear. Note that it is not possible to erase the user program if the protection code for the user program is enabled in the System Menu.
- <LOG>** The ProStar 430 keeps a log of system-relevant events (**<EVENTS>**); records error messages that have been generated) and keeps count of actions of valves and syringe movements (**<COUNT>**). A message appears after every 50,000 syringe actions and after every 200,000 syringe valve actions: "Lifetime of syringe (valve) maybe exceeded. Check for possible leakage!"
- Syringe: if you do not replace the syringe at this moment and tell the system "not to display this message again", the message will not be displayed again until 50,000 more syringe actions have been counted. Syringe valve: will have to be replaced and the counter for valve actions will have to be reset by the maintenance engineer.
- <DEFAULT ALL>** to change all software settings to default. All series, methods, templates and the user program (unless protected by protection code) will be erased.

NOTE: If **<DEFAULT ALL>** is selected, check whether hardware configuration still is compatible with settings entered in the System Menu.

- <SSV> (option)** Use this key to start a procedure in which all lines of the solvent selection valve can be primed. The menu offers the following possibilities:
- <SSV1>** to open the corresponding port of the solvent selection valve.
to
- <SSV6>**
- <PRIME>** to prime the selected solvent line with one syringe volume. The last selected port will remain active after leaving the SSV mode.

<COOL> (option) Use this key to enter the programming mode for Peltier plate cooling. The programmable temperature range is 4°C to 40°C. The maximum cooling capacity is approximately 12°C below ambient (refer to the Appendix for specifications). Connect the condensed water and leakage connector to a waste container on the floor to drain condensed water.

If the cool option is switched **<ON>** the following soft function keys can be selected:

<MANUAL> temperature control will remain OFF until it is switched on again by the user (in this menu).

<AUTOMATIC> Temperature control will be switched OFF after all programmed series have been executed.

<DATE-TIME> Temperature control will be switched OFF at a date and time that can be programmed.

<SERIAL> Use this key to put the ProStar 430 in serial mode to allow for control of the autosampler by way of PC (RS232 interface). Select a device identifier in the System Menu. If a method protection code was defined in the system settings, this code must be entered to get access to serial mode. The following soft function keys appear:

<PANIC> Press this key to begin a stop sequence in which all tubing is rinsed and the valve and I/O ports are reset. At the end of the sequence serial mode is resumed.

<EXIT> Press this key to end serial mode and return to the Ready Menu.

<SERVICE> Indicates service required. To be used by authorized personnel only. The Service Menu is protected by a service code.

System Menu

The System Menu contains the following soft function keys:

<GENERAL> Press this key to enter values for:

Loop Volume, Needle Tubing, Syringe Volume	These values have to be entered because the ProStar 430 can be fitted with various types of syringes. Every needle and needle tubing volume requires a different minimum flush volume. The default flush volume equals two times the volume of needle and tubing.
---	---

NOTE:	The 10 mL syringe can only be selected if the use of the User Program is enabled in the System Menu (Usage Menu). When the 10 mL syringe is selected, only the User Program can be used. Use of other types of methods will be disabled.
--------------	--

Syringe Speed and Scale Factor	The aspirating speed of the syringe used in injection methods can be adapted depending on viscosity of samples. Alternatively syringe speeds can be reduced by entering a scale factor. The syringe speed will be the scale factor multiplied by the syringe speed. The speed of the syringe during the wash or the rinsing procedure of the buffer are not affected by this setting.
Needle Height	Distance between the needle point and the plate holder (default: 5 mm). The value is only used in injection methods (for mix methods this value is programmable in the method itself).
Skip Missing Vials	Appears only if a plate with 48 vials is selected in System Menu (Plates Menu). YES means that empty spaces are skipped during the run. NO means that the ProStar 430 will stop if an empty space is encountered during the run; an error code will be generated.
Air Segment	Decide whether an air segment will be used for analytical runs.
Headspace Pressure	To switch headspace pressure on or off. The ProStar 430 uses headspace pressure to facilitate transport of sample into the loop. The compressor will always be used during a wash procedure. Please note that accuracy and reproducibility may decrease if headspace pressure is switched off. However, headspace pressure will only be useful if sample wells are airtight.
Time Display	Offers a choice between two types of time representation.
Key click, error beep and alarm buzzer	Offer the possibility to switch sound signals on or off.

Table 7 Overview of General Settings

General setting for	Default	Possible range
loop volume	100 µL	5 - 1000 µL
needle tubing	15 µL	1 - 200 µL
syringe volume	250 µL	100, 250, 500, 1000 µL or 10 mL
syringe speed	normal	low, normal or high
scale factor	1.0	0.1 - 1.0
needle height	5 mm	0 - 40 mm
skip missing vials	yes	yes or no
air segment	yes	yes or no
headspace pressure	no	yes or no
time base display	HH:MM:SS	H:MM:SS or H:MM:mm
key click	on	on or off
error beep	on	on or off
alarm buzzer	on	on or off

<USAGE>

Press this key to enter the following Usage settings:

	SETTINGS	DEFAULT
Protection Code	To enter a code for protection of methods. Enter a six digit code (000000-999999) for protection of all methods. Press CL to erase the code. If a code has been defined it is not possible to enter the System Menu and the programming menus without entering the protection code.	None
Timebase methods	To enable or disable the possibility to program timebase methods. The ProStar 430 controls other connected equipment during analysis time. Program timebase methods in the Methods Menu	Disabled

	SETTINGS	DEFAULT
Mix methods	To enable or disable the possibility to program mix methods for the ProStar 430. Program mix methods in the Methods Menu. Please note that the ProStar 430 cannot analyze priority samples during a run if the mix method is enabled.	Disabled
User program	To enable or disable the possibility to program a user program. If this function is enabled it is possible to enter a user program protection code (6 digits). Program the user program in the Methods Menu. Please note that the ProStar 430 cannot analyze priority samples during a run if the user program is enabled.	Disabled
Labeled wells	To enable or disable the possibility to program labeled wells. Program the location of labeled wells in the Series Menu.	Disabled
Templates	To enable or disable the possibility to program templates. Program templates in the Methods Menu.	Disabled
Calibration well	To enable or disable the possibility of programming calibration wells. Program the location of the calibration wells in the Series Menu.	Disabled

NOTE: You are advised to disable as many functions in the Usage Menu as possible to make sure that other menus do not contain possibilities that are irrelevant for the type of analyses you are to perform.

<PLATES>	Use this key to define the type of plates that will be used. Five types can be selected: 96-low (default), 96-high, 384-low, 48-vials, or user defined. After a plate type has been selected, enter:
Well processing method	In rows (left to right) or in columns (top to bottom).
First transport vial	Enter a number 1 – 4, or press CL.

	<p>Last transport vial</p> <p>Enter a number 1 - 4.</p> <p>Vials can be placed in any of the four positions. Transport vials must be placed in a continuous row.</p>
<IO>	<p>Use this key to enter the I/O configuration mode and define the following:</p>
	<p>Inject-marker pulse length</p> <p>To define the length of the inject-marker pulse.</p>
	<p>Well-marker pulse length</p> <p>To define the length of the well-marker pulse.</p>
	<p>Labeled well marker pulse length</p> <p>To define the length of the well-marker pulse of the labeled well.</p>
	<p>Input edge next injection</p> <p>To define the edge sensitive inputs for the next injection.</p>
	<p>Input edge next well</p> <p>To define the edge sensitive inputs for the next well.</p>
	<p>Freeze input active</p> <p>To define whether the freeze input is active when high, or freeze input is active when low.</p>
	<p>Reset outputs after last series</p> <p>To determine whether the outputs should be reset to default after the last series.</p> <p>Refer to page 35 for more specific information on I/O connections.</p>

Table 8 Overview I/O Menu

Possibility	Default	Range
inject-marker pulse length	1.0 s	0.1 - 2.0
well-marker pulse length	1.0 s	0.1 - 2.0
labeled well marker pulse length	1.0 s	0.1 - 2.0
input edge next injection	falling	falling or rising
input edge next well	falling	falling or rising
freeze input active	low	low or high
reset outputs after last series	no	yes or no

<CLOCK>

Use this key to switch the system clock on or off. Select ON to enter the clock menu in which you can set date (yy,mm,dd) and time (hh,mm). This date and time will be displayed in the Ready Menu.

<COMM.>

Use this key to define a device identifier for communication with other equipment (e.g., a PC). An identifier between 20 and 29 can be selected for the ProStar 430.

Methods Menu

This menu allows the user to program various types of methods: it is possible to define 24 separate injection methods, 5 wash methods, 5 timebase methods, 9 mix methods and one user program.

It is also possible to program a combination of methods and save them in a template. The settings entered in the System Menu determine the possibilities offered by the Methods Menu.

<TEMPLATE>

Use this key to enter a menu in which the contents of a template can be defined. First assign a number to the template, then link the numbers of methods to the template. A maximum of 24 templates can be programmed. The following items can be entered to fill a template:

User program instead of methods

If soft function key **<YES>** is selected, the complete template is filled with the user program; no other methods can be added. If soft function key **<NO>** is selected the template can be filled with the following:

- **mix method number**
- **injection method number**
- **wash method number**
- **timebase method number**

<METHODS>

Use this key to enter a menu in which methods can be defined:

<MIX>

(If enabled in System Menu, Usage Menu)

Use this key to program a method that allows you to perform pre-injection sample handling, e.g. pre-column derivatization, dilution or adding of internal standard. Nine mix methods can be programmed; the maximum number of steps that can be programmed for the total of 9 mix methods and the user program is 240. Assign a number to the mix method. The Mix Menu appears:

<EDIT>: an existing step or a new step for a new mix method

<INSERT>: a new step in an existing method before the displayed step

<DELETE>: the displayed step

"End of mix method" means that the mix method is empty; if an existing mix method is selected, the first line of the mix method is displayed. Scroll through the steps of the existing method with the cursor keys and use the soft function keys to enter changes in an existing method.

The following types of steps can be programmed for a mix method:

- <ASPIRATE>** (Sample, air, destination, reagent A-D) a programmed volume. Speed of syringe can be selected from 1-9. Height (H) indicated is the distance of the needle point to the plate holder (default: 5 mm). The maximum amount which can be aspirated is the total volume of the syringe.
- <DISPENSE>** (Sample, waste, destination, reagent A-D) a programmed volume from the buffer tubing. Speed of the syringe can be selected from 1 - 9. Height (H) indicated is the distance of the needle point to the plateholder (default: 5 mm). It is possible to dispense a larger volume than the volume aspirated in previous actions. The aspirated amount will be complemented with liquid from the wash solvent bottle to total the programmed dispense volume.
- <WAIT>** To define a pause (H:MM:SS, maximum of 9 hours, 59 minutes and 59 seconds).
- During the pause, the needle will move to home position (if the previous step is an aspirate or dispense action). If you want the needle to stay in the same position, an aspirate or dispense step of 0 μ L must be programmed at the desired position
- <REPEAT>** Enter the number of steps that must be repeated and how often they must be repeated.
- <WASH>** Enter the volume for needle wash. Buffer is rinsed to waste.
- <INJECTION>** Use this key to program a method that defines injection methods (max. 24) for a run of the types full loop, partial loopfill or μ L pick-up. Enter a number for the injection method you are going to program. If the selected method is locked because of changes in the settings (System Menu) after programming the method, the word LOCK is displayed. The method can be unlocked by programming valid values in the method itself or by restoring the values in the System Menu.
- Use the soft function keys to select an injection method, then enter values for:

Analysis time	The time between switching the injection valve to inject and the start of processing the next sample.
Flush volume	The amount of sample taken from a vial before the loop is filled with sample. Default value: 30 µL (combined with an air segment).

NOTE: Flush volumes of less than twice the volume of the needle and tubing will result in decreased performance.

Number of injections per vial	Maximum value is 9.
Injection volume	Can be entered for each injection per vial. The maximum programmable injection volumes are:
<i>Partial loopfill</i>	50% of the programmed loop volume
<i>µL Pick-up</i>	50% of the programmed loop volume Not programmable, is equal to the loop volume but needs more sample to fill the loop (3 x loop volume for loop volumes < 100 µL; 2 x loop volume for loop volumes ≥ 100 µL - 499 µL; 1.5 x loop volume for loop volumes ≥ 500 µL).
<i>Full Loop</i>	Not programmable, is equal to the loop volume but needs more sample to fill the loop (3 x loop volume for loop volumes < 100 µL; 2 x loop volume for loop volumes ≥ 100 µL - 499 µL; 1.5 x loop volume for loop volumes ≥ 500 µL).

<WASH> Use this key to program wash methods. It is possible to program a wash between injections, samples or series. For each wash method the volume of wash solvent can be defined. The minimum programmable volume is 300 µL.

<TIMEBASE> (If enabled in System Menu, Usage Menu)

Press this key to enable control of the optional ISS valve and other devices via auxiliary or binary outputs. A maximum of 5 timebase methods can be programmed. The menu offers the following soft function keys:

<AUX> Scroll through all program lines by pressing **E** or select AUX to move to the next auxiliary.

<VALVES>

Controls the ISS valve and the solvent selection valve. The ISS valve can only be programmed if the optional ISS valve is installed. 6-1 and 2-1 refer to the interconnected ports of the valves. Press **E** to scroll through programming lines (Only if SSV option is installed). Enter the time and the SSV port number (value between 1 and 6).


Figure 11 ISS Valve, interconnected ports

<CODE>

Outputs: enter a time and a value between 1 and 15, hexadecimal output. Press **E** to scroll through the programming lines.

<END>

To enter the end time for timed events program; press **E** to scroll through the programming lines. If no value is filled in or if **CL** is pressed, the ProStar 430 will automatically generate an end time. The end time is equal to the analysis time programmed in the injection method used in the same series.

If end time exceeds the programmed analysis time, this end time overrules the analysis time.

It is possible to program events after the end time, but these events are not carried out during a run.

<USER PROGRAM>

The user program offers the possibility to program all possible actions required for a sample handling sequence in separate steps. Note that the total number of steps for the user program and all nine mix methods cannot exceed 240. The user program can be protected by a special user program protection code (System Settings, Usage Menu). If no user program has been programmed yet, "end of user program" is displayed. Otherwise, the first line of the programmed method appears. The soft function keys appear:

<EDIT>: an existing step or program a new step for the user program.

<INSERT>: a new step before the displayed step

<DELETE>: delete the displayed step

The edit and insert menus offer the following soft function keys:

<ASPIRATE> a programmed volume from sample well, ambient air, destination vial, wash, or one of the reagent vials into the buffer tubing. Speed and height of syringe can be entered (refer to Table 9). The maximum volume that can be aspirated is the total volume of the syringe.

<DISPENSE> a programmed volume from the buffer tubing into the sample well, waste, destination vial, wash or one of the reagent vials. Speed and height of syringe can be entered (refer Table 9).

It is not possible to dispense a larger volume than the total volume aspirated in previous actions.

<SYR VALVE> to control the connections of the syringe to one of its three tubes:

<NEEDLE>: connection to sample needle

<WASH>: connection to wash solvent bottle

<WASTE>: connection to syringe waste tubing.

<SYR> to control the movements of the syringe.

<LOAD>: the syringe with the programmed volume

<UNLOAD>: the syringe with the programmed volume

<HOME>: the volume previously aspirated will be dispensed to the last programmed position, and the syringe will be initialized again.

<WASH> to execute a needle wash; the content of the buffer tubing is not rinsed to waste before the start of the wash. The programmed volume of wash solvent is used to wash the needle at the wash position.

The wash position may be contaminated with the contents of the buffer tubing, which may generate cross-contamination. To prevent contamination of the wash position, program a dispense to waste action before programming a wash action.

- <VALVES>** To program positions of high pressure valves (ISS, injector valve, SSV). The injector valve has two positions: **<INJECT>** and **<LOAD>**. The ISS optional valve has positions 1-6 and 1-2 (Figure 11).
- <WAIT>** to program a pause (max. 9 hours, 59 minutes, 59 seconds).

NOTE: During the pause, the needle will move to home position (if the previous step is an aspirate or dispense action). If you want the needle to stay in the same position, an aspirate or dispense step of 0 µL must be programmed at the desired position.

- <COMPRES>** to activate the compressor to put air pressure on a sample. The compressor will stay active until it is switched off (in a next programmed step). The compressor will be automatically switched off at the end of the needle wash routine if a needle wash is used.
- <AUX>** to control the four standard auxiliaries (contact closures).
- <WAIT-IN>** to program a pause in which the ProStar 430 waits for one of the four inputs to become **<HIGH>** or **<LOW>** before continuing with the next step.
- <PROG-OUT>** to define two programmable outputs (contact closures). These are similar to the auxiliaries, but only available in the user program.
- <CODE>** to program the output to the connector P3 TIMED OUTPUTS. This is a HEX output in the range 0 to 15.
- <MARKERS>** The markers normally generated in the ProStar 430 are not active in the user program, but can be programmed in this screen. Select marker and status (inject, vial, labeled).
- <SSV>** (option): to define the Solvent Selection Valve (SSV) port position, range 1 to 6.

Table 9 Syringe Speed

SPEED	SYRINGE				
	100 μ L	250 μ L	500 μ L	1000 μ L	10 mL
1	50 μ L/min	125 μ L/min	250 μ L/min	500 μ L/min	0.5 mL/min
2 (low)	125 μ L/min	315 μ L/min	630 μ L/min	1255 μ L/min	1.3 mL/min
3 (normal)	250 μ L/min	625 μ L/min	1250 μ L/min	2495 μ L/min	2.5 mL/min
4 (high)	375 μ L/min	940 μ L/min	1880 μ L/min	3765 μ L/min	3.8 mL/min
5	770 μ L/min	1920 μ L/min	3840 μ L/min	7680 μ L/min	7.7 mL/min
6	1070 μ L/min	2675 μ L/min	5335 μ L/min	10665 μ L/min	10.7 mL/min
7	1370 μ L/min	3430 μ L/min	6855 μ L/min	13710 μ L/min	13.7 mL/min
8	1745 μ L/min	4365 μ L/min	8725 μ L/min	17450 μ L/min	17.5 mL/min
9	2135 μ L/min	5335 μ L/min	10670 μ L/min	21335 μ L/min	21.3 mL/min

During the dispense action the pressure in the buffer tubing will increase. To prevent damage of the buffer tubing, the flow should not exceed the value of 6 mL/min for water. (Maximum speed 9 for 100 μ L and 250 μ L-syringes, speed 6 for a 500 μ L syringe and speed 4 for 1000 μ L and 10 mL syringe.) If more viscous liquids are used the speeds should be reduced.

Series Menu

This menu allows you to define the run sequence in a series. A maximum of 24 series can be programmed. A series contains information about the methods to be used for a range of wells. This can be a template, a separate method (mix, injection, wash, timebase), or the user program. Information on location of wells, labeled wells or calibration wells is also programmed in a series.

NOTE:	The settings entered in the System Menu and the methods defined in the Methods Menu determine which possibilities appear in the Series Menu.
-------	--

Table 10 Series Parameters

<i>Without templates</i>	<i>With templates</i>
<input type="radio"/> Use user program Yes/No <input checked="" type="radio"/> Injection method number <input checked="" type="radio"/> Wash method number <input type="radio"/> Time base methods number <input type="radio"/> Mix method number <i>Time base and mix method are only available if enabled in the System Menu</i>	<input checked="" type="radio"/> Template number
<input type="radio"/> Use calibration wells Yes/No <ul style="list-style-type: none"> <input type="radio"/> First calibration well <input type="radio"/> Last calibration well <input type="radio"/> No. of samples between calibration <i>Calibration wells are only available if enabled in the System Menu; not available if Mix Method has been programmed</i>	
	<input checked="" type="radio"/> First sample well <input checked="" type="radio"/> Last sample well
	<i>Only if a mix method has been programmed, or if used in user program:</i> <ul style="list-style-type: none"> <input type="radio"/> First destination well <input type="radio"/> Vial Reagent-A <input type="radio"/> Vial Reagent-B <input type="radio"/> Vial Reagent-C <input type="radio"/> Vial Reagent-D
	<i>Only if the use of labelled wells has been enabled in the System Menu:</i> <ul style="list-style-type: none"> <input type="radio"/> Labelled well no. 1 <input type="radio"/> Labelled well no. 2 <input type="radio"/> Labelled well no. 3 <input type="radio"/> Labelled well no. 4

☒ marked questions are always asked in series

☐ marked questions depend on the used methods and the settings entered in the System Menu.

Explanation

After you have entered the required settings in the System Menu and after you have programmed methods to be used for an analytical run, you can press Series to enter the Series Menu. Table 10 gives an overview of the items you have to define for the Series.

With Templates

If you are going to execute an analytical run by way of a template, you will only be asked to enter the template number and to indicate the location of the first sample well and the last sample well.

Without Templates

If you are going to execute an analytical run without using a template, you will be asked to enter an injection method number and a wash method number, and you will have to indicate the location of the first sample well and last sample well.

If you have enabled use of calibration wells in the System Menu (Usage Menu), you will have to define whether you will use calibration wells, and indicate the location of the first and last calibration well, and indicate the number of wells between calibration wells.

However, if you have for example enabled use of a Mix Method in the System Menu (Usage Menu), you will also have to define the location of the First destination well and Reagent vials.

NOTE: Series are stored in the ProStar 430 memory for as long as the power is on. As soon as power is switched off, all programmed series will be deleted.

It is not possible to leave the Series Menu before all values have been programmed.


Figure 12 Injection sequence with 3 calibration wells between every 5 wells


Programming Options

[illegible]

Programming Options

	Ready Menu						◀ [Menu] ▶					
	<Exchange>	<Wash>	<Sry End>	<ULTS>		<SERIAL>	<SSV>**	<COOL>**	<SERVICE>			
				<Copy>	<Erase>	<Log>	<Default All>					
	Move Plates Move Lifts Adjust Plate Height	Perform standard wash	Move syringe	Copy method	Erase methods & Userprog	Log-book	All methods erased, all settings to default	Comm. protocol	Prime solvent valve	Cool ON/OFF Setpoint Switch ON/OFF		
										Service Mode		
[System]	System Settings										◀ [Menu] ▶	
	<GENERAL>		<USAGE>		<PLATES>		<IO>		<CLOCK>		<COMM.>	
	Loop volume Needle tubing Syringe volume Needle Height Skip missing wells Air segment Headspace pressure Time base display Key click Error buzzer Alarm buzzer		Method protection code Use time base method Use mix method Use user-program User-program protection code Use labeled wells Use templates Use calibration wells		Type of plate Process in ROWS/ COLUMNS Number of transport vials		Inj. Marker pulse Well marker pulse Lab. Well marker pulse Input edge next inj. Input edge next vial Freeze input active Reset outputs after series		System clock Date Time		Device identifier	
Method	METHODS PROGRAMMING											
	<TEMPLATE>		<METHOD>							<USERPROG>		
	User prog* or Mix method injection Wash method Timebase method*		METHODS <INJECTION>				<WASH>		<TIMEBASE>*		<MIX>*	
			Injection Methods				Wash volume Wash between: –Series –Wells –Injections		Timebase program: 4x AUX 1 4x AUX 2 4x AUX 3 4x AUX 4 4x ISS 8x SSV 8x CODE OUT END TIME		Mix method: ASPIRATE –sample –destination –air –reagent DISPENSE –sample –destination –waste –reagent REPEAT Wait Wash	
			<FULL>				<PARTIAL>		<PICK-UP>			
			Analysis time Flush volume Injection/well				Analysis time Flush volume Injection/well Inj. volumes		Analysis time Injection/well Inj. volumes			
[Series]	With Templates Template Use calibration wells (yes/no) * – calibration wells – calibration interval Sample Wells Destination wells * Reagent wells * Labeled wells *		Without Templates User program (yes/no) or Mix method* Injection method Wash method Timebase method* Use calibration wells (yes/no)* – calibration wells – calibration interval Sample Wells Destination wells * Reagent wells * Labeled wells *									
[Start/ Stop]	First series Last series <START> <REMOTE>											
[Priority]	Template Well number		Injection Method Wash method Timebase method Well Number									
[Interrupt]	Not Used											
[Help]	Context and status sensitive help screen (if available)											

Maintenance and Troubleshooting

Analytical problems

In case of analytical problems you will have to determine whether they are caused by the autosampler or by the rest of the system.

NOTE:

Quick check!

Replace the valve with a manual injection valve to discriminate between valve problems and other problems.

Perform a number of manual Full loop injections. If the results are good, the autosampler is faulty; if not, the HPLC system should be checked.




Figure 13 Fluidics Diagram

In the flowcharts on the next pages it is assumed that the ProStar 430 does not display any error messages.

Please keep in mind that analytical problems might be caused by external influences, like temperature and/or light-sensitive samples. For this reason it is important to make sure that the application was running without problems before and that no changes have been made in the settings (System Menu).

Error Codes

The ProStar 430 will display an **error message** if the user tries to enter invalid data. Information on the allowed range will be displayed. If something goes wrong in the physical operations of the ProStar 430, an **error code** will be displayed. Press the **Start/Stop** key twice to lift the message and try to repair the failure condition with the help of the explanation of the code concerned.

Injection Valve and ISS Unit

ERROR 11	Injection valve is not in a valid position.
ERROR 12	The injection valve did not switch within 1.5 seconds.
ERROR 13	The switching time of the injection valve exceeds 500 msec.
ERROR 14	ISS valve is not in a valid position.
ERROR 15	The ISS valve did not switch within 1.5 seconds.

Syringe Dispenser Unit

ERROR 21	The syringe valve did not switch.
ERROR 22	The syringe did not reach home position in time.
ERROR 23	The syringe spindle did not make the correct number of rotations.
ERROR 24	The spindle does not rotate.
ERROR 25	The syringe valve did not find a valid position.

Injection Needle Unit

ERROR 30	The sample needle arm did not reach or leave home position (vertical).
ERROR 31	The sample needle arm is in an invalid horizontal position while moving down.
ERROR 32	The sample needle arm did not reach or leave destination within a certain time (horizontal).
ERROR 34	Sample needle arm not in vertical home position while moving horizontally.

ERROR 39	Vial sensor sticks.
ERROR 40	The sample needle spindle does not rotate correctly.
ERROR 41	The sample needle did not reach or leave home position.
ERROR 42	The sample needle is not at home position.
ERROR 53	The sample needle arm is not in the home position while moving the plate.

Plate

ERROR 59	Missing plate.
----------	----------------

Vials

ERROR 60	Missing vial. Only available when Skip Missing Vial is set to NO in the System Settings and during the execution of the Mix of a sample on 48-vial plate.
ERROR 62	Missing transport vial.
ERROR 64	Missing vial for reagent A.
ERROR 65	Missing vial for reagent B.
ERROR 66	Missing vial for reagent C.
ERROR 67	Missing vial for reagent D.
ERROR 68	Missing destination vial.
ERROR 69	Not enough transport liquid available due to missing transport vials.

Electronics

ERROR 71	Flexprint of the sample needle is not connected.
ERROR 72	Invalid configuration of the ProStar 430, PCB missing.
ERROR 73	Current limit of the external I/O exceeded.
ERROR 75	Error occurred during initialization, the ProStar 430 cannot start.

Plate holder

ERROR 90	Plate home time-out, plate did not reach home position (home error).
ERROR 91	Plate did not reach or leave home position during run.
ERROR 92	Plate holder missing.
ERROR 93	Dirt on the plate holder.

Appendix

Specifications

General

Working temperature	5 - 40°C
Storage temperature	-25 / +60°C
Humidity	20 - 80% RH
Sample viscosity	0.1 - 5 cP
Sample capacity	1 plate + 4 x 10 mL vials.
Plates (refer to figure B.1)	types: 48-vials, 96-low, 96-high, 384-low Max. height: 47 mm (incl. cap); min. height: 12.5 mm (incl. cap). Note: Select high (height per plate 42 mm) or low (height per plate 14.5 mm) plate in System Menu.
Vial dimensions	6 x 8 vials, dimensions base plate 128 x 85,8 mm Maximum outer diameter for each of 48-vials: 12 mm
Reagent/transport vials	Dimensions 10 mL vials: 22 x 47 mm (incl. cap); min. 40 mm (incl. cap)
Loop volume	5-1000 µL
Puncturing accuracy needle	± 0.5 mm
Pre-puncturing septa/caps	With air needle, dual needle action
Plate & height detection	Missing plate & height detection by vial sensor

Plate tray drive	> 5 cm/sec
Switching time injection valve	Electrically < 100 msec
Headspace pressure	Built-in compressor
Wash solvent	250 mL internal wash solvent bottle
Dispenser syringe	100 µL , 250 µL (standard), 500 µL, 1000 µL, or 10 mL
Wetted parts	SS316, PTFE, TEFZEL, VESPEL, Glass, Teflon (Optional PEEK, Fused silica)

Analytical performance

Capped and Seal Vials

Reproducibility

full loop injections	RSD \leq 0.3 %
partial loopfill injections	RSD \leq 0.5 %, injection volumes > 5 µL, with headspace pressure on the vial and 30 µL pre-flush with air segment.

Open vials/plates, valid at 1.0 cP:

Reproducibility

full loop injections	RSD \leq 0.3 %
partial loopfill injections	RSD \leq 0.5 %; injection volumes > 5 µL, without headspace pressure on the vial and 30 µL pre-flush with air segment.
µL pick-up injections	RSD \leq 1.0 %; injection volumes > 5 µL. without headspace pressure on the vial.

Memory Effect

< 0.01 % with programmable needle wash

Programming

Injection methods	Full loop injections Partial loopfill injections
Injection volume	5 µL - 1 mL, full loop, depending on system settings 1 - 500 µL, with 0.1 µL increments for partial loopfill 1 µL - max. volume, with 0.1 µL increments for µL pick-up max. volume = ((loop volume - 3 x needle volume)/2)
Injections per vial	Max. 9 (volumes are programmable for each injection)
Analysis time	Max. 9 hrs. 59 min. 59 sec.
Needle wash	Programmable (between injections, wells or series)
Priority sample	Freely programmable
Series	Freely programmable, 24 series max.

Physical

Dimensions (W x D x H)	280 mm x 400 mm x 440 mm
Weight	23 kg 26 kg if cooling option is installed

Electrical

Power requirements	115 Vac; + 15/-20 %; 50 Hz/60 Hz; 250 VA 230 Vac; + 15/-20 %; 50 Hz/60 Hz; 250 VA
Fuses	For 115 Vac; two 5.0 T-Type fuses (¼" x 1¼", UL/CSA) For 230 Vac; two 2.5 AT-fuses (5 x 20 mm, IEC 127) All fuses UL-listed and CSA-certified

Communication

OUTPUTS

Inject marker
Well marker
Labeled well marker
Stop I/O
4 Auxiliary outputs
2 Programmable outputs
Alarm output
4 Bit timebase

INPUTS

Next injection input
Next well input
Freeze input
Stop I/O
4 Programmable inputs

Options

Sample cooling

Built-in Peltier cooling processing unit
Programmable Range : 4° C - 40° C
Cooling capacity: Ambient -12° C (measured on cooling plate)

Integrated stream switching (ISS)

One optional high pressure 6-port valve

Solvent selection valve (SVV)

6-way integrated solvent selection valve

Options and Accessories

The following accessories are available for the ProStar 430:

Part Number	Description
R000088050	Needle wash insert vial
R000088051	Needle body
0391478800	Operation Manual
0392611060	Valco Injection Valve
0392611061	Valco Biocompatible Injection Valve
0392607967	Flanged tube fitting 1/8 (pck/5)
0392607968	Flanged tube fitting 1/16 (pck/5)
0392607831	Plunger replacement tip 100 μ L (pck/10)
0392607832	Plunger replacement tip 500 μ L (pck/10)
0392607964	Plunger replacement kit 1 mL (pck/10)
R000088052	Sample needle including tubing ferrules and nuts
R000088064	Serum sample needle including tubing ferrules and nuts
R000088068	Peek sample needle including tubing ferrules and nuts
0392607835	Silica-coated/steel sample needle including tubing ferrules and nuts
0392607836	Silica-coated/steel serum needle including tubing ferrules and nuts
R000088065	Sample needle extra long including tubing ferrules and nuts
0392607837	Serum sample needle extra long including tubing ferrules and nuts
R000088066	Peek sample needle extra long including tubing ferrules and nuts
0392607838	Fused silica needle including tubing ferrules and nuts
0392607839	Fused silica needle extra long including tubing ferrules and nuts
0392607804	Buffer tubing 500 μ L

Part Number	Description
0392607805	Buffer tubing 2000 µL
0392607971	Syringe waste tubing
0392607972	Syringe wash tubing
0392607847	Buffer tubing 15 mL
R000088053	Air/prepuncturing needle
0392607963	Plunger replacement tip 250 µL (pck/10)
0392607848	Shipkit ProStar 430
0392607908	Syringe valve
0392607922	Luerlock connection syringe valve
R000088055	Wash solvent bottle 250 mL
0392607962	Syringe 1000 µL
R000088067	Syringe 10 mL
R000088057	Syringe 100 µL
0392607810	Syringe 250 µL
R000088059	Syringe 500 µL
0190018200	5 µL SS Loop for Valco Valve
0190018202	20 µL SS Loop for Valco Valve
0392611011	50 µL SS Loop for Valco Valve
0190018204	100 µL SS Loop for Valco Valve
0392611012	250 µL SS Loop for Valco Valve
0190018205	500 µL SS Loop for Valco Valve
0392611013	1000 µL SS Loop for Valco Valve
0392611014	2000 µL SS Loop for Valco Valve
0392611015	5000 µL SS Loop for Valco Valve
0392611016	10 mL SS Loop for Valco Valve
0392611017	20 µL PEEK Loop for Valco Valve

Part Number	Description
0392611018	50 µL PEEK Loop for Valco Valve
0190018208	100 µL PEEK Loop for Valco Valve
0392611019	250 µL PEEK Loop for Valco Valve
0392611020	500 µL PEEK Loop for Valco Valve
0392611021	1000 µL PEEK Loop for Valco Valve
0392611022	2000 µL PEEK Loop for Valco Valve
2869450201	Package of 10 Valco Ferrules
2869450101	Package of 10 Valco Nuts
0392607983	Valco Fingertight PEEK Nut
0392607984	Valco Fingertight PEEK Ferrule

The following Plates and Inserts are available for the ProStar 430.

Part Number	Description
0392607760	350 µL Round Well Titer Plate Polypropylene 50/box
0392607761	1 mL Round Well Titer Plate Polypropylene, 50/box
0392607762	350 µL Square Well Titer Plate Polypropylene, 50/box
0392607763	1 mL Square Well Titer Plate Polypropylene, 50/box
0392607764	2 mL Square Well Titer Plate Polypropylene, 50/box
0392607765	Piercible Round Well Sealing Cap Polypropylene, 50/box
0392607766	Piercible Square Well Sealing Cap Polypropylene, 50/box
0392607767	300 µL Glass Insert with Snap Cap for 03-926077-61, 100/box
0392607768	300 µL Glass Insert for 03-926077-61, 100/box
0392607769	Snap Caps for Glass Inserts, 100/box
0392607770	384 Well Polypropylene Plates 50/Box