

7 14.48 - 14.7 - 15.1


MSD Tune Report, G1946: MM-APCI Positive Mode - Standard Scan

Instrument : Instrument 1


Rev. B.02.01-SR1 [260]

Sat Jan 07 06:07:10 2012

C:\Chem32\1\1946DTUN\ATUNES.TUN



Scan: 50.00 - 2300.00 Samples: 12 Thresh: 0 Step: 0.10
22501 peaks Base: 922.10 Abundance: 472064 F Raw,Gauss 0.3



m/z	gf	Abund	Rel Abund	Iso m/z	Iso Abund	Iso Ratio
121.06		56080	100.00	122.04	4473	7.98
622.04		314560	560.91	623.10	45808	14.56
922.10		472064	841.77	923.10	98176	20.80
1521.95		169536	302.31	1523.00	58264	34.37
2121.90		28048	50.01	2122.90	12555	44.76

```


Mode:MM-APCI   Polarity: POS



Fragment      VAR WidthGain  820
Skim1         20 WidthOffs   VAR
Lens1         2.0 MassGain  74.15
Lens2DC       20.0 MassOffs   VAR
Iris          350
               Energy       5.0
               OpolPeak     300
               OpolKnee     312
Gain          1.00 QuadDC    0.00


EMV           1400 Samples   12
VCap          2500 Averages   3
Corona         6.0 StepSize  0.10

StepSize      0.10
CapCur        31 DryingGas  5.0
ChamCur       2.41 Gas Temp  350
CoronaVol     4961 Vaporizer  201
Quad Temp      99 Neb Pres   20
RoughVac       2.89 HighVac  1.6e-5


```







Scan: 50.00 - 2300.00 Samples: 12 Thresh: 0 Step: 0.10
 22501 peaks Base: 806.00 Abundance: 273152 F Raw,Gauss 0.3



m/z	gf	Abund	Rel Abund	Iso m/z	Iso Abund	Iso Ratio
119.05		22496	100.00	120.02	1756	7.81
555.96		161472	717.78	556.94	18776	11.63
806.00		273152	1214.22	807.00	47736	17.48
1306.00		163968	728.88	1307.00	48328	29.47
1806.00		27640	122.87	1806.90	8792	31.81

MSD Tune Report, G1946: MM-APCI Negative mode - Standard Scan

Instrument : Instrument 1

Rev. B.02.01-SR1 [260]

Sat Jan 07 06:12:12 2012


C:\Chem32\1\1946DTUN\ATUNES.TUN

```
Mode:MM-APCI   Polarity: NEG

Fragment  VAR WidthGain  821
Skim1     20 WidthOffs  VAR
Lens1    -2.4 MassGain  74.25
Lens2DC   24.5 MassOffs  0.668
Iris      -25
          Energy      5.0
          OpolPeak    250
          OpolKnee     280
Gain      1.00 QuadDC   0.00

EMV       1400 Samples   12
VCap      2500 Averages   1
Corona     6 StepSize  0.10

StepSize  0.10
CapCur    16 DryingGas  5.0
ChamCur   2.09 Gas Temp  350
CoronaVol  3494 Vaporizer 200
Quad Temp  99 Neb Pres   20
RoughVac   2.88 HighVac  1.6e-5
```



MSD Tune Report, G1946: MM-APCI Negative mode - Standard Scan

Instrument : Instrument 1

Rev. B.02.01-SR1 [260]

Sat Jan 07 06:12:12 2012

C:\Chem32\1\1946DTUN\ATUNES.TUN

